

# ***Service Quick Reference***

---

---

## **HANDBOOK**

---

---

**JUNE 1, 2005**

CONTENTS SUBJECT TO CHANGE WITHOUT NOTICE

# **STERIS®**



STERIS Corporation • STERIS Engineering Service  
2424 West 23rd Street • Erie, PA 16506

# ◆ INDEX ◆

	<u>PAGE</u>
• SERVICE KIT LISTING	1
• P.M. PACKAGES/RE-NEW PACKS	
Polaris Lights	5
2080M Tables	6
• P.M. PACKS	
Steam Sterilizers & Steam Generators	7
Low Temperature Sterilization &	
EO Products	22
Decontamination	24
Surgical Tables	27
Surgical Lights	27
Specialty Surgical Products	28
Patient Handling & Transport Systems	28
Miscellaneous Equipment	29
• LAMPS	30
• STERILIZABLE/DISPOSABLE SURGICAL LIGHT HANDLES	31
• 2000 SERIES STERILIZERS W/PRINTCON P/C BOARDS	32
• PRINTED CIRCUIT BOARDS	33
• EAGLE 3000 STAGE 3 CONTROLS	36
• PC BOARDS RETURN LIST	37
• TABLE LIFT CYLINDERS/ ACCESSORIES	38
• 3080/3085 HYDRAULIC CYLINDERS	39
• HAUSTED CYLINDER ASSEMBLIES, COLUMNS & MEDIJACKS	40
• SURGICAL AND ACCESSORY PADS	41

**NOTE:** Parts that are marked **NLA** indicate that the part is **No Longer Available**.

# INDEX (Cont'd.)

	<u>PAGE</u>
• BRISTOL TO U.E. RECORDER/ CONTROLLER REPLACEMENT CHART	45
• REMANUFACTURED STERILIZER DOORS	47
• DOOR GASKETS	48
• VALVE CHART FOR EAGLE 2000 STAGE 2 STERILIZERS	50
• VALVE CHART FOR EAGLE 3000 STAGE 2 STERILIZERS	52
• VALVE CHART FOR RENAISSANCE STERILIZERS	54
• VALVE CHART FOR EAGLE 3000 STAGE 3 STERILIZERS	55
• VALVE CHART FOR 2400 & 3400 AMSCOMATIC WASHER/STERILIZER	57
• SOLENOID VALVES FOR SYSTEM 2S/LAB 90, CENTURY/MILLENNIUM 116/120, AND AMSCO CENTURY MEDIUM STERILIZERS	59
• SMALL & MEDIUM STERILIZER SAFETY VALVES	63
• SAFETY VALVE SPECIFICATIONS	64
• EO EAGLE MONITOR S/S FITTINGS & COMPONENTS	65
• SURGICAL PRODUCTS CONSUMABLES	66
• TOOL LIST	70
• MISCELLANEOUS TOOLS	71
• EO/SAFETY RELATED TOOLS	71

## INDEX (Cont'd.)

	<u>PAGE</u>
• SERVICE LITERATURE	72
• STERIS TECHNICIAN LABELS, TAGS & SIGNS	72
• PRINTER PAPER MATRIX	73
• SMALL & MEDIUM STERILIZERS VACUUM BREAKERS/REPAIR KITS	74
• ELECTRO-STEAM STEAM GENERATOR HEATING ELEMENTS	75
• CHROMALOX STEAM GENERATOR HEATING ELEMENTS	78
• SERVICE SUPPORT RESOURCES - ERIE CONTACTS	80
• KEY ACCOUNTS	81
• NOTES	83

## KIT LISTING Sterilizers

Part Number	Description
P056401-052 .....	Kit, Deletion of Manual Exhaust, All Eagle Sterilizers
P129363-893 .....	Kit, Data Acq., VHP 1000
P129363-910 .....	Kit, Humidity Sensor, VHP 1000
P129363-958 .....	Kit, Impact Printer, VHP 1000
P150822-426 .....	Kit, Air Inlet Valve Replacement, Eagle 48" with Envirogard
P150822-442 .....	Kit, Air Inlet Valve Replacement, Eagle Small with Envirogard
P150822-443 .....	Kit, Air Inlet Valve Replacement, Eagle Small Printcon
P150822-492 .....	Kit, Air Inlet Valve Replacement, Eagle Small w/o Envirogard
P150822-568 .....	Kit, Air Inlet Valve Replacement, Eagle Sterilizer Aerator
P150822-569 .....	Kit, Air Inlet Valve Replacement, Eagle Sterilizer Aerator
P150822-570 .....	Kit, Air Inlet Valve Replacement, Eagle Sterilizer Aerator
P150828-322 .....	Kit, Drain Line Restrictor, Medallion Eagle Sterilizers
P760266-023 .....	Kit, Auto Flush & Drain, Chromalox Stand-alone Steam Generators - CH05 thru CH14
P761951-001 .....	Kit, Liquid Sensing S.D. Medallion W/S
P762650-131* .....	Kit, Steam to Electric, 20" Medallion GP/Lab/Lab ISO (20kW,208/240Volt,3PH)
P762650-132* .....	Kit, Steam to Electric, 20" Medallion GP/Lab/Lab ISO (20kW,440/480Volt,3PH)
P762650-136* .....	Kit, Steam to Electric, 16" Medallion GP/Lab/Lab ISO (20kW,208/240Volt,3PH)
P762650-137* .....	Kit, Steam to Electric, 16" Medallion GP/Lab/Lab ISO (20kW,440/480Volt,3PH)
P762650-147* .....	Kit, Steam to Electric, 16" Medallion GP/Lab/Lab ISO (20kW,208/240Volt,1PH)
P762650-150* .....	Kit, Steam to Electric, 20" Medallion GP/Lab/Lab ISO (20kW,380/415Volt,3PH)
P762650-151* .....	Kit, Steam to Electric, 20" Medallion GP/Lab/Lab ISO (20kW,575/600Volt,3PH)
P762650-152* .....	Kit, Steam to Electric, 20" Medallion GP/Lab/Lab ISO (20kW,208/240Volt,1PH)
P762650-154* .....	Kit, Steam to Electric, 16" Medallion GP/Lab/Lab ISO (20kW,380/415Volt,3PH)
P762650-155* .....	Kit, Steam to Electric, 16" Medallion GP/Lab/Lab ISO (20kW,575/600Volt,3PH)
P762650-368 .....	Kit, Steam to Electric, 16" & 20" Eagle Century Grav/Prevac/ISO (30kw,380/415Volt,3PH)
P763856-091 .....	Kit, Insulating & Fastening - 36x48x60 S.D./ DBL
P764184-091 .....	Kit, Insulation for Small Medallion/Eagle
P764315-666 .....	Kit, Drain Pan - Medallion PIWS
P764320-235* .....	Kit, Steam to Electric, 20" Eagle 2000 (20kW,208/240Volt)
P764320-236* .....	Kit, Steam to Electric, 20" Eagle 2000 Grav/Grav ISO (20kW,440/480Volt,3PH)
P764320-603* .....	Kit, Steam to Electric, 16" Eagle 2000 (20kW,440/480Volt)
P764320-668 .....	Kit, Water Line Conversion - Small Eagle GP
P764320-669 .....	Kit, Water Line Conversion - Med Eagle GP
P764320-670 .....	Kit, Water Line Conversion - Small Eagle VAC
P764320-671 .....	Kit, Water Line Conversion - Med Eagle VAC
P764320-782* .....	Kit, Steam-to-Electric Conversion - 16 Eagle 2000 (20kW,208/240Volt)
P764321-115 .....	Kit, Pre-Heat Bypass

P764322-890*	Kit, Steam to Electric, 20" Eagle 2000 Grav/Grav ISO (20kW,380/415Volt,3PH)
P764323-140	Kit, Steam to Electric, 20" Stage 2 Grav (Flex Conduit) (30kW,208/240Volt)
P764323-224*	Kit, Steam to Electric, 16" Eagle 2000 Grav/Grav ISO (20kW,575/600Volt,3PH)
P764323-225*	Kit, Steam to Electric, 20" Eagle 2000 Grav/Grav ISO (20kW,575/600Volt,3PH)
P764323-963	Kit, Capacitor Replacement, Stage 2
P764323-985	Kit, Steam to Electric, 20" E 3000 Stage 2, Grav w/Anti-Flood (Neoprene Cables) (30kW,208/240Volt,3PH)
P764323-986	Kit, Steam to Electric, 20" E 3000 Stage 2, Grav w/Anti-Flood (Neoprene Cables) (30kW,440/480Volt,3PH)
P764323-987	Kit, Steam to Electric,16" E 3000 Stage 2 Grav (Neoprene Cables) (30kW,208/240Volt,3PH)
P764323-988	Kit, Steam to Electric,16" E 3000 Stage 2 Grav (Neoprene Cables) (30kW,440/480Volt,3PH)
P764324-014	Kit, Amscomatic 2400 Float Type Water Level Sensor
P764324-085	Kit, Temp Cntrl Replace, 3 Temp Switch, 20x20 Cryo (78378-091)
P764324-089	Kit, Steam to Electric,16" E 3000 Stage 2 Grav (Flex Conduit) (30kW,208/240Volt,3PH)
P764324-099	Kit, Steam to Electric,16" E 3000 Stage 2 Grav (Flex Conduit) (30kW,440/480Volt,3PH)
P764324-100	Kit, Steam to Electric, 20" E 3000 Stage 2, Grav w/Anti-Flood (Flex Conduit) (30kW,440/480Volt,3PH)
P764324-101	Kit, Steam to Electric,16" E 3000 Stage 2 Vac w/Anti-Flood (Neoprene Cables) (30kW,208/240Volt,3PH)
P764324-102	Kit, Steam to Electric,16" E 3000 Stage 2 Vac w/Anti-Flood (Neoprene Cables) (30kW,440/480Volt,3PH)
P764324-103	Kit, Steam to Electric, 20" E 3000 Stage 2, Vac w/Anti-Flood (Neoprene Cables) (30kW,208/240Volt,3PH)
P764324-104	Kit, Steam to Electric, 20" E 3000 Stage 2, Vac w/Anti-Flood (Neoprene Cables) (30kW,440/480Volt,3PH)
P764324-105	Kit, Steam to Electric,16" E 3000 Stage 2 Grav (Neoprene Cables) (30kW,575/600Volt,3PH)
P764324-106	Kit, Steam to Electric,16" E 3000 Stage 2 Grav (Flex Conduit) (30kW,575/600Volt,3PH)
P764324-107	Kit, Steam to Electric, 20" E 3000 Stage 2, Grav w/Anti-Flood (Neoprene Cables) (30kW,575/600Volt,3PH)
P764324-108	Kit, Steam to Electric, 20" E 3000 Stage 2, Grav w/Anti-Flood (Flex Conduit) (30kW,575/600Volt,3PH)
P764324-109	Kit, Steam to Electric,16" E 3000 Stage 2 Vac w/Anti-Flood (Neoprene Cables) (30kW,575/600Volt,3PH)
P764324-110	Kit, Steam to Electric,16" E 3000 Stage 2 Vac w/Anti-Flood (Flex Conduit) (30kW,575/600Volt,3PH)
P764324-111	Kit, Steam to Electric, 20" E 3000 Stage 2, Vac w/Anti-Flood (Neoprene Cables) (30kW,575/600Volt,3PH)
P764324-112	Kit, Steam to Electric, 20" E 3000 Stage 2, Vac w/Anti-Flood (Flex Conduit) (30kW,575/600Volt,3PH)
P764324-115	Kit, Steam to Electric,16" E 3000 Stage 2 Vac w/Anti-Flood (Flex Conduit) (30kW,208/240Volt,3PH)
P764324-116	Kit, Steam to Electric,16" E 3000 Stage 2 Vac w/Anti-Flood (Flex Conduit) (30kW,440/480Volt,3PH)
P764324-117	Kit, Steam to Electric, 20" E 3000 Stage 2, Vac w/Anti-Flood (Flex Conduit) (30kW,208/240Volt,3PH)

P764324-118 .....	Kit, Steam to Electric, 20" E 3000 Stage 2, Vac w/Anti-Flood (Flex Conduit) (30kW,440/480Volt,3PH)
P764324-121 .....	Kit, Steam to Electric,16" E 3000 Stage 2 Grav (Flex Conduit) (30kW,380/415Volt,3PH)
P764324-122 .....	Kit, Steam to Electric,16" E 3000 Stage 2 Grav (Neoprene Cables) (30kW,380/415Volt,3PH)
P764324-123 .....	Kit, Steam to Electric,16" E 3000 Stage 2 Vac w/Anti-Flood (Flex Conduit) (30kW,380/415Volt,3PH)
P764324-124 .....	Kit, Steam to Electric,16" E 3000 Stage 2 Vac w/Anti-Flood (Neoprene Cables) (30kW,380/415Volt,3PH)
P764324-125 .....	Kit, Steam to Electric, 20" E 3000 Stage 2, Grav w/Anti-Flood (Flex Conduit) (30kW,380/415Volt,3PH)
P764324-126 .....	Kit, Steam to Electric, 20" E 3000 Stage 2, Grav w/Anti-Flood (Neoprene Cables) (30kW,380/415Volt,3PH)
P764324-127 .....	Kit, Steam to Electric, 20" E 3000 Stage 2, Vac w/Anti-Flood (Flex Conduit) (30kW,380/415Volt,3PH)
P764324-128 .....	Kit, Steam to Electric, 20" E 3000 Stage 2, Vac w/Anti-Flood (Neoprene Cables) (30kW,380/415Volt,3PH)
P764324-572 .....	2400/3400 Trap Conversion Parts List
P764326-277 .....	Kit, Steam to Electric, 20" E 3000 Stage 3 Grav/Vac (30kW,380/415Volt,3PH)
P764326-287 .....	Kit, Snap Switch Model 53 Cyclomatic/Medallion X Replacement Kit
P764326-466 .....	Kit, STM to ELC Eagle 3000 Stage 3 16x16" (208/240V, 3-110V,1)
P764326-467 .....	Kit, STM to ELC Eagle 3000 Stage 3 16x16" (440/480V, 3-110V,1)
P764326-468 .....	Kit, STM to ELC Eagle 3000 Stage 3 20x20" (208/240V, 3-110V,1)
P764326-469 .....	Kit, STM to ELC Eagle 3000 Stage 3 20x20" (440/480V, 3-110V,1)
P764326-719 .....	Kit, Steam to Electric, 16" E 3000 Stage 3 Grav/Vac (30kW,380/415Volt,3PH)
P764326-876* .....	Kit, Steam to Electric, 20" Eagle 2000 Grav/Grav ISO (30kW,380/415Volt,3PH)
P764327-571 .....	Kit, Steam to Electric, 16" Eagle Century Grav/Prevac/ISO (30kW,240Volt,3PH)
P764327-572 .....	Kit, Steam to Electric, 20" Eagle Century Grav/Prevac/ISO (30kW,240Volt,3PH)
P764327-573 .....	Kit, Steam to Electric, 16" Eagle Century Grav/Prevac/ISO (30kW,440/480Volt,3PH)
P764327-574 .....	Kit, Steam to Electric, 20" Eagle Century Grav/Prevac/ISO (30kW,440/480Volt,3PH)
P764328-028 .....	Kit, Auto Flush & Drain, Integral Steam Generators
P764328-068 .....	Kit, Lovato Contactor, Steam Generator
P764327-571 .....	Kit, Steam to Electric, 16" Century (30kW,240Volt,3PH)
P764327-572 .....	Kit, Steam to Electric, 20" Century (30kW,240Volt,3PH)
P764327-573 .....	Kit, Steam to Electric, 16" Century (30kW,440/480Volt,3PH)
P764327-574 .....	Kit, Steam to Electric, 20" Century (30kW,440/480Volt,3PH)
P762650-368 .....	Kit, Steam to Electric, 16", 20" Century (30kW,380/415Volt,3PH)
P764328-397 .....	Kit, Steam to Electric, 16" Eagle Century Grav/Prevac/ISO (30kW,208Volt,3PH)
P764328-398 .....	Kit, Steam to Electric, 20" Eagle Century Grav/Prevac/ISO (30kW,208Volt,3PH)

- P764328-425\* ..... Kit, Steam to Electric, 20" Eagle 2000 Grav/Grav ISO  
(30kW,208/240Volt,3PH)
- P764328-494\* ..... Kit, Steam to Electric, 20" Eagle 2000 Grav/Grav ISO  
(30kW,440/480Volt,3PH)

*\* Order Water Sensing Kit if kit (anti-flood) has been previously added to Medallion or Eagle 2000 Series units.*

### **Tables/Lights**

- P129360-660 ..... Kit, 3080 Table, Andrews Frame Kit
- P129360-764 ..... Kit, Hand Control - 2080RC
- P150823-309 ..... Kit, Back Section Rail Extension, 3080
- P630939-039 ..... Kit, AC to Battery Conversion 3080
- P758132-091 ..... Kit, Low Leakage Light - 9' Single Track
- P764320-733 ..... Kit, Conversion, 2080M to 2080MIA
- P764321-139 ..... Kit, Conversion, 2080 L&RC to 2080 LIA&RCIA
- P764323-356 ..... Kit, Conversion - 9' Split Track Conv Polaris/Challenge
- P764328-571 ..... Kit, Cord Replacement-Hand Control, 3080/3085 RC/RL/SP
- P764326-652 ..... Kit, Examiner 10 to Single Ceiling Mount
- P764326-653 ..... Kit, Examiner 10(2) to Dual Ceiling Mount
- P764326-654 ..... Kit, Examiner 10 Light to 54, Open Track Mount
- P764326-655 ..... Kit, Examiner 10 Light to 108", Open Track
- P764326-689 ..... Kit, 3080/3085 Foot Control Cord Replacement
- P764326-690 ..... Kit, 3080/3085 Foot Control Switch Replacement
- P764326-694 ..... Kit, Retrofit Plunger, 62" I.A. Board

### **Allied Products**

- B602047-467 ..... Kit, Upgrade, Thermostat, 4th Only Elec
- P150823-295 ..... Kit, Soap pumps, Uni/ Flexmatic
- P150825-321 ..... Kit, Service-Sonic Heater Conversion
- P754372-001 ..... Kit, Whistle Valve, 3/8" - Aeroflush
- P762327-001 ..... Kit, Second Heater - Amsco 18, Warming Cabinet, 4th Gen.
- P762328-001 ..... Kit, Second Heater - Amsco 24" Warm Cabinet/Lower 4th Gen
- P762329-001 ..... Kit, Heater - 18" Warm Cab, Lower Compartment 4th Gen
- P762330-001 ..... Kit, Heater - 24" Warm Cab, Lower Compartment 4th Gen
- P764315-573 ..... Kit, Soap Pump Replacement Unimatic Only
- P764315-805 ..... Kit, Conversion - Adj Thermostat, Lower Compartment 4th Gen
- P764316-785 ..... Kit, Thermostat - Medallion Warming Cabinet 3rd Gen
- P764317-910 ..... Kit, Convert Steam to Electric Warming Cabinet/Upper  
Compartment 4th Gen
- P764321-303 ..... Kit, Thermometer (all Warming Cabinets)
- P764323-919 ..... Kit, Scrub Sink, Field Install (Double) Uniflex
- P764326-570 ..... Kit, Timed Injection, Basil
- P764326-571 ..... Kit, Conductivity Injection, Basil
- P764326-602 ..... Kit, Scrub Sink Timer, Flexmatic



**RE-NEW PACK  
POLARIS LIGHT CENTRA MOUNT  
PARABOLIC REFLECTOR A**

**Re-New Pack - Polaris Light Centramount Parabolic Reflector A -  
P764325-813**

Rack, Support .....	P092918-002
Retro-Reflector Assembly, Scorp .....	P093900-819
Cap Assembly .....	P141211-042
Spacer .....	P150543-002
Expansion Plug .....	P150621-001
Spring Washer .....	P150622-001 (2)
Screw, 6-32 Hex Soc Fl Hd .....	P150824-132 (10)
Kit, Brush Holder Assembly .....	P764315-406
Kit, Lens, Polaris Light .....	P764317-855
Kit, Roundel Filter & Parabolic - Pair "A" .....	P056938-805
Screw .....	P150544-001
Gasket, Filter .....	P150824-297

**RE-NEW PACK  
POLARIS LIGHT CENTRA MOUNT  
PARABOLIC REFLECTOR B**

**Re-New Pack Polaris Light Centramount Parabolic Reflector B -  
P764325-814**

Rack, Support .....	P092918-002
Retro-Reflector Assembly, Scorp .....	P093900-819
Cap Assembly .....	P141211-042
Spacer .....	P150543-002
Expansion Plug .....	P150621-001
Spring Washer .....	P150622-001 (2)
Screw, 6-32 Hex Soc Fl Hd .....	P150824-132 (10)
Kit, Brush Holder Assembly .....	P764315-406
Kit, Lens, Polaris Light .....	P764317-855
Kit, Roundel Filter & Parabolic - Pair "B" .....	P056938-806
Screw .....	P150544-001
Gasket, Filter .....	P150824-297

# **RE-NEW PACK** **2080 M TABLE AFTER 8/74**

## **Re-New Pack 2080 M Table after 8/74 - P764325-822**

Taper Pin .....	P015342-045
Spring .....	P047777-091
Shift Lever Assembly .....	P055093-091
Caster Assembly .....	P056397-135 (4)
X-Ray Top Spacer .....	P077643-056 (2)
Drive Shaft .....	P078222-061
Socket Head Screw .....	P080001-091
Spring .....	P080117-091
Bearing .....	P082633-001 (4)
Trend Crank .....	P092013-001
Drive Crank .....	P092016-001
X-Ray Top Spacer .....	P092475-001 (2)
Kidney Handle Assembly .....	P092588-001
Screw .....	P150055-001 (2)
Pin and Foot Assembly .....	P150202-001 (4)

## **Re-New Pack 2080 M Table Lift Cylinder/Pump Replacement after 8/74 - P764325-823**

Pump Assembly .....	P055087-001
Lift Cylinder .....	P136375-001
Hydraulic Oil .....	P764322-636

## **P.M. Pack - Miscellaneous**

*TLT Pad Set .....	P150830-177
*Armboard Repair Kit .....	P764316-574

*\*Order Individual Parts as Required*

## P.M. PACKS

### STEAM STERILIZERS AND STEAM GENERATORS

<u>Description</u>	<u>Part Number</u>
<b>Eagle 1016 &amp; Eagle 10 Tabletop Sterilizers (PMCL P764327-012)</b>	
Major Insp. #1 - order PM Pack	P764330-224
<b>FLJ Vacumatic Floor Loading Sterilizer w/Stage 3 Control (PMCL P764327-013)</b>	
Major Insp. #1 - order PM Pack	P764330-268
<b>FLS Vacumatic Floor Loading Sterilizer w/Stage 3 Control (PMCL P764327-013)</b>	
Major Insp. #1 - order PM Pack	P764330-269
<b>Chromalox CES CH05 through CH08 - Carbon Steel only (30 kW through 75 kW) (PMCL P764327-014)</b>	
Major Insp. #1 - order Descaler* 105108 and PM Pack	P764330-355
Major Insp. #2 - order PM Pack	P764330-356
<b>Chromalox CES CH10 through CH14 - Carbon Steel only (100 kW through 180 kW) (PMCL P764327-014)</b>	
Major Insp. #1 - order Descaler* 105108 and PM Pack	P764330-363
Major Insp. #2 - order PM Pack	P764330-364
<b>AMSCO Century Medium 26x37.5 Prevac &amp; SFPP Sterilizer (PMCL P764327-265)</b>	
Major Insp. #1	PM Pack not required
Major Insp. #2 - order PM Pack	P764332-098
<b>Century 16"/20" Gravity &amp; Vacumatic Steam Sterilizers - Single Door, Steam (PMCL P764327-263)</b>	
Major Insp. #1	PM Pack not required
Major Insp. #2 - order PM Packs	P764326-475 (1) P764326-480 (1) P764326-619 (1) P764332-065 (1)
<b>Century 16"/20" Gravity &amp; Vacumatic Steam Sterilizers - Single Door, Electric (PMCL P764327-263 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #2 for Generator - order PM Pack	P764330-247
Major Insp. #2 for Sterilizer - order PM Packs	P764326-475 (1) P764326-480 (1) P764326-619 (1) P764332-065 (1)

\*Descaler 105108 (consumable)

## **P.M. PACKS (continued)**

### **STEAM STERILIZERS AND STEAM GENERATORS (continued)**

<u>Description</u>	<u>Part Number</u>
<b>Century 16"/20" Gravity &amp; Vacumatic Steam Sterilizers - Double Door, Steam (PMCL P764327-263)</b>	
Major Insp. #1	PM Pack not required
Major Insp. #2 - order PM Packs	P764326-475(1) P764326-480(1) P764326-619(1) P764332-065(1)
<b>Century 16"/20" Gravity &amp; Vacumatic Steam Sterilizers - Double Door, Electric (PMCL P764327-263 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer	PM Pack not required
Major Insp. #2 for Generator - order PM Pack	P764330-247
Major Insp. #2 for Sterilizer - order PM Packs	P764326-475(1) P764326-480(1) P764326-619(1) P764332-065(1)
<b>Millennium 16"/20" Gravity &amp; Vacumatic Steam Sterilizers - Single Door, Steam (PMCL P764327-263)</b>	
Major Insp. #1	PM Pack not required
Major Insp. #2 - order PM Packs	P764326-475(1) P764326-480(1) P764326-619(1) P764332-065(1)
<b>Millennium 16"/20" Gravity &amp; Vacumatic Steam Sterilizers - Single Door, Electric (PMCL P764327-263 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer	PM Pack not required
Major Insp. #2 for Generator - order PM Pack	P764330-247
Major Insp. #2 for Sterilizer - order PM Packs	P764326-475(1) P764326-480(1) P764326-619(1) P764332-065(1)
<b>Millennium 16"/20" Gravity &amp; Vacumatic Steam Sterilizers - Double Door, Steam (PMCL P764327-263)</b>	
Major Insp. #1	PM Pack not required
Major Insp. #2 - order PM Packs	P764326-475(1) P764326-480(1) P764326-619(1) P764332-065(1)

*\*Descaler 105108 (consumable)*

## P.M. PACKS (continued)

### STEAM STERILIZERS AND STEAM GENERATORS (continued)

<u>Description</u>	<u>Part Number</u>
<b>Millennium 16"/20" Gravity &amp; Vacumatic Steam Sterilizers - Double Door, Electric (PMCL P764327-263 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer	PM Pack not required
Major Insp. #2 for Generator - order PM Pack	P764330-247
Major Insp. #2 for Sterilizer - order PM Packs	P764326-475 (1) P764326-480 (1) P764326-619 (1) P764332-065 (1)
<b>Eagle 180 Steam Sterilizer - Steam (PMCL P764327-092)</b>	
Major Insp. #1 - order following part(s) separately	
• Gasket, Door 12" sq.	P136800-905
• Grease, 4 oz Tube w/Nozzle	P764322-635
Major Insp. #2 - order PM Pack	P764330-277
<b>Eagle 180 Steam Sterilizer - Electric (PMCL P764327-092)</b>	
Major Insp. #1 - order following part(s) separately	
• Gasket, Door 12" sq.	P136800-905
• Grease, 4 oz Tube w/Nozzle	P764322-635
Major Insp. #2 - order PM Pack	P764330-278
<b>System 2S Steam Sterilizer (PMCL P764327-120)</b>	
Major Insp. #1	PM pack not required
Major Insp. #2 - order following part(s) separately	
• Kit, Valve Repair 1/4" - S1	P764328-740
• Kit, Valve Repair 1/4" - S2, S35, S37	P764328-732 (3)
• Kit, Valve Repair 1/2" - S3	P764328-738
• Kit, Valve Repair 1/2" - S7, S50	P764328-736 (2)
• Kit, Repair 1/2" - TR1	P762326-001
• Valve, Swg Ck 1/4" - CK1, CK2, CK4	P129376-130 (3)
• Valve, Check 3/8" NPT - CK3	P005424-091
• Valve, Swing Check 1/4" NPT - CK5	P010278-091
• Valve, Ck SS Swaglock 1/4" NPT-CK6	P338521-266
• Filter, Element Air - AF1	P101006-172
• Valve, Safety 1/2" 50 psig - RV1	P093918-013
<b>Eagle 2000 16" Vacumatic Steam Sterilizer - Steam (PMCL P764327-259)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-002
Major Insp. #2 - order PM Pack	P764332-062

\*Descaler 105108 (consumable)

**P.M. PACKS (continued)**

**STEAM STERILIZERS AND STEAM GENERATORS (continued)**

<u>Description</u>	<u>Part Number</u>
--------------------	--------------------

**Eagle 2000 20" Vacumatic Steam Sterilizer - Steam  
(PMCL P764327-259)**

Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 - order PM Pack	P764332-062

**Eagle 3000 Stage 2 16" Vacumatic Steam Sterilizer - Steam  
(PMCL P764327-259)**

Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-002
Major Insp. #2 - order PM Pack	P764332-062

**Eagle 3000 Stage 2 20" Vacumatic Steam Sterilizer - Steam  
(PMCL P764327-259)**

Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 - order PM Pack	P764332-062

**Eagle 3000 Stage 3 16" Vacumatic Steam Sterilizer - Steam  
(PMCL P764327-259)**

Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-002
Major Insp. #2 - order PM Pack	P764332-058

**Eagle 3000 Stage 3 20" Vacumatic Steam Sterilizer - Steam  
(PMCL P764327-259)**

Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 - order PM Pack	P764332-058

**Renaissance 16" Vacumatic Steam Sterilizer - Steam  
(PMCL P764327-259)**

Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-002
Major Insp. #2 - order PM Pack	P764332-058

**Renaissance 20" Vacumatic Steam Sterilizer - Steam  
(PMCL P764327-259)**

Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 - order PM Pack	P764332-058

## P.M. PACKS (continued)

### STEAM STERILIZERS AND STEAM GENERATORS (continued)

<u>Description</u>	<u>Part Number</u>
<b>Eagle 2000 16" Vacumatic Steam Sterilizer - Single Door, Electric w/Chromalox (CHS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package	P754360-002
Major Insp. #2 for Generator - order PM Pack	P764330-247
Major Insp. #2 for Sterilizer - order PM Pack	P764332-062
<b>Eagle 2000 16" Vacumatic Steam Sterilizer - Single Door, Electric w/Electrosteam Solid State (AS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package	P754360-002
Major Insp. #2 for Generator - order PM Pack	P764330-250
Major Insp. #2 for Sterilizer - order PM Pack	P764332-062
<b>Eagle 2000 16" Vacumatic Steam Sterilizer - Single Door, Electric w/Electrosteam Non-Solid State (AS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-246
Major Insp. #1 for Sterilizer - order Parts Package	P754360-002
Major Insp. #2 for Generator - order PM Pack	P764330-251
Major Insp. #2 for Sterilizer - order PM Pack	P764332-062
<b>Eagle 2000 20" Vacumatic Steam Sterilizer - Electric w/Chromalox (CHS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 for Generator - order PM Pack	P764330-247
Major Insp. #2 for Sterilizer - order PM Pack	P764332-062

*\*Descaler 105108 (consumable)*

## P.M. PACKS (continued)

### STEAM STERILIZERS AND STEAM GENERATORS (continued)

<u>Description</u>	<u>Part Number</u>
<b>Eagle 2000 20" Vacumatic Steam Sterilizer - Electric w/Electrosteam Solid State (AS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 for Generator - order PM Pack	P764330-250
Major Insp. #2 for Sterilizer - order PM Pack	P764332-062
<b>Eagle 2000 20" Vacumatic Steam Sterilizer - Electric w/Electrosteam Non-Solid State (AS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-246
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 for Generator - order PM Pack	P764330-251
Major Insp. #2 for Sterilizer - order PM Pack	P764332-062
<b>Eagle 3000 Stage 2 16" Vacumatic Steam Sterilizer - Single Door, Electric w/Chromalox (CHS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package	P754360-002
Major Insp. #2 for Generator - order PM Pack	P764330-248
Major Insp. #2 for Sterilizer - order PM Pack	P764332-062
<b>Eagle 3000 Stage 2 16" Vacumatic Steam Sterilizer - Single Door, Electric w/Electrosteam (AS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package	P754360-002
Major Insp. #2 for Generator - order PM Pack	P764330-249
Major Insp. #2 for Sterilizer - order PM Pack	P764332-062
<i>*Descaler 105108 (consumable)</i>	



**P.M. PACKS (continued)**

**STEAM STERILIZERS AND STEAM GENERATORS (continued)**

<u>Description</u>	<u>Part Number</u>
<b>Eagle 3000 Stage 2 20" Vacumatic Steam Sterilizer - Electric w/Chromalox (CHS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 for Generator - order PM Pack	P764330-248
Major Insp. #2 for Sterilizer - order PM Pack	P764332-062
<b>Eagle 3000 Stage 2 20" Vacumatic Steam Sterilizer - Electric w/ Electrosteam (AS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 for Generator - order PM Pack	P764330-249
Major Insp. #2 for Sterilizer - order PM Pack	P764332-062
<b>Eagle 3000 Stage 3 16" Vacumatic Steam Sterilizer - Single Door, Electric w/Chromalox (CHS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package	P754360-002
Major Insp. #2 for Generator - order PM Pack	P764330-247
Major Insp. #2 for Sterilizer - order PM Pack	P764332-058

*\*Descaler 105108 (consumable)*

## P.M. PACKS (continued)

### STEAM STERILIZERS AND STEAM GENERATORS (continued)

<u>Description</u>	<u>Part Number</u>
<b>Eagle 3000 Stage 3 20" Vacumatic Steam Sterilizer - Electric w/ Chromalox (CHS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 for Generator - order PM Pack	P764330-247
Major Insp. #2 for Sterilizer - order PM Pack	P764332-058
<b>Renaissance 16" Vacumatic Steam Sterilizer - Single Door, Electric w/ Chromalox (CHS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package	P754360-002
Major Insp. #2 for Generator - order PM Pack	P764330-247
Major Insp. #2 for Sterilizer - order PM Pack	P764332-058
<b>Renaissance 20" Vacumatic Steam Sterilizer - Electric w/Chromalox (CHS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 for Generator - order PM Pack	P764330-247
Major Insp. #2 for Sterilizer - order PM Pack	P764332-058
<b>Eagle 2000 16" Gravity Steam Sterilizer - Steam (PMCL P764327-259)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-002
Major Insp. #2 - order PM Pack	P764332-061
<b>Eagle 2000 20" Gravity Steam Sterilizer - Steam (PMCL P764327-259)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 - order PM Pack	P764332-061

\*Descaler 105108 (consumable)

## P.M. PACKS (continued)

### STEAM STERILIZERS AND STEAM GENERATORS (continued)

<u>Description</u>	<u>Part Number</u>
<b>Eagle 3000 Stage 2 16" Gravity Steam Sterilizer - Steam (PMCL P764327-259)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-002
Major Insp. #2 - order PM Pack	P764332-061
<b>Eagle 3000 Stage 2 20" Gravity Steam Sterilizer - Steam (PMCL P764327-259)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 - order PM Pack	P764332-061
<b>Eagle 3000 Stage 3/Renaissance 16" Gravity Steam Sterilizer - Steam (PMCL P764327-259)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-002
Major Insp. #2 - order PM Pack	P764332-057
<b>Eagle 3000 Stage 3/Renaissance 20" Gravity Steam Sterilizer - Steam (PMCL P764327-259)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 - order PM Pack	P764332-057
<b>Eagle 2000 16" Gravity Steam Sterilizer - Electric w/Chromalox (CHS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-002
Major Insp. #2 for Generator - order PM Pack	P764330-247
Major Insp. #2 for Sterilizer - order PM Pack	P764332-061
<b>Eagle 2000 16" Gravity Steam Sterilizer - Electric w/Electrosteam Solid State (AS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-002
Major Insp. #2 for Generator - order PM Pack	P764330-250
Major Insp. #2 for Sterilizer - order PM Pack	P764332-061

\*Descaler 105108 (consumable)

## P.M. PACKS (continued)

### STEAM STERILIZERS AND STEAM GENERATORS (continued)

<u>Description</u>	<u>Part Number</u>
<b>Eagle 2000 16" Gravity Steam Sterilizer - Electric w/Electrosteam Non-Solid State (AS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-246
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-002
Major Insp. #2 for Generator - order PM Pack	P764330-251
Major Insp. #2 for Sterilizer - order PM Pack	P764332-061
<b>Eagle 2000 20" Gravity Steam Sterilizer - Electric w/Chromalox (CHS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 for Generator - order PM Pack	P764330-247
Major Insp. #2 for Sterilizer - order PM Pack	P764332-061
<b>Eagle 2000 20" Gravity Steam Sterilizer - Electric w/Electrosteam Solid State (AS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 for Generator - order PM Pack	P764330-250
Major Insp. #2 for Sterilizer - order PM Pack	P764332-061

*\*Descaler 105108 (consumable)*

## P.M. PACKS (continued)

### STEAM STERILIZERS AND STEAM GENERATORS (continued)

<u>Description</u>	<u>Part Number</u>
<b>Eagle 2000 20" Gravity Steam Sterilizer - Electric w/Electrosteam Non-Solid State (AS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-246
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 for Generator - order PM Pack	P764330-251
Major Insp. #2 for Sterilizer - order PM Pack	P764332-061
<b>Eagle 3000 Stage 2 16" Gravity Steam Sterilizer - Single Door, Electric w/Chromalox (CHS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package	P754360-002
Major Insp. #2 for Generator - order PM Pack	P764330-248
Major Insp. #2 for Sterilizer - order PM Pack	P764332-061
<b>Eagle 3000 Stage 2 16" Gravity Steam Sterilizer - Single Door, Electric w/Electrosteam (AS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package	P754360-002
Major Insp. #2 for Generator - order PM Pack	P764330-249
Major Insp. #2 for Sterilizer - order PM Pack	P764332-061

*\*Descaler 105108 (consumable)*

## P.M. PACKS (continued)

### STEAM STERILIZERS AND STEAM GENERATORS (continued)

<u>Description</u>	<u>Part Number</u>
<b>Eagle 3000 Stage 2 20" Gravity Steam Sterilizer - Electric w/Chromalox (CHS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 for Generator - order PM Pack	P764330-248
Major Insp. #2 for Sterilizer - order PM Pack	P764332-061
<b>Eagle 3000 Stage 2 20" Gravity Steam Sterilizer - Electric w/Electrosteam (AS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 for Generator - order PM Pack	P764330-249
Major Insp. #2 for Sterilizer - order PM Pack	P764332-061
<b>Eagle 3000 Stage 3 16" Gravity Steam Sterilizer - Single Door, Electric w/Chromalox (CHS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package	P754360-002
Major Insp. #2 for Generator - order PM Pack	P764330-247
Major Insp. #2 for Sterilizer - order PM Pack	P764332-057

*\*Descaler 105108 (consumable)*

## P.M. PACKS (continued)

### STEAM STERILIZERS AND STEAM GENERATORS (continued)

<u>Description</u>	<u>Part Number</u>
<b>Eagle 3000 Stage 3/Renaissance 20" Gravity Steam Sterilizer - Electric w/ Chromalox (CHS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 for Generator - order PM Pack	P764330-247
Major Insp. #2 for Sterilizer - order PM Pack	P764332-057
<b>Renaissance 16" Gravity Steam Sterilizer - Single Door, Electric w/ Chromalox (CHS) (PMCL P764327-259 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package	P754360-002
Major Insp. #2 for Generator - order PM Pack	P764330-247
Major Insp. #2 for Sterilizer - order PM Pack	P764332-057
<b>Eagle 2000 24x36" Vacumatic Steam Sterilizer (PMCL P764327-261)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P764324-196
Major Insp. #2 - order PM Pack	P764332-060
<b>Eagle 3000 Stage 2 24x36" Vacumatic Steam Sterilizer (PMCL P764327-261)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P764324-196
Major Insp. #2 - order PM Pack	P764332-060
<b>Eagle 3000 Stage 3/Renaissance 24x36" Vacumatic Steam Sterilizer (PMCL P764327-261)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P764324-196
Major Insp. #2 - order PM Pack	P764332-056
<b>Eagle 2000 24x36" Gravity Steam Sterilizer (PMCL P764327-261)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P764324-196
Major Insp. #2 - order PM Pack	P764332-059

*\*Descaler 105108 (consumable)*

## P.M. PACKS (continued)

### STEAM STERILIZERS AND STEAM GENERATORS (continued)

<u>Description</u>	<u>Part Number</u>
<b>Eagle 3000 Stage 2 24x36" Gravity Steam Sterilizer (PMCL P764327-261)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P764324-196
Major Insp. #2 - order PM Pack	P764332-059
<b>Eagle 3000 Stage 3/Renaissance 24x36" Gravity Steam Sterilizer (PMCL P764327-261)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P764324-196
Major Insp. #2 - order PM Pack	P764332-055
<b>Medallion 20x20" Vacumatic S Steam Sterilizer (PMCL P764327-155)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 - order PM Pack	P764330-273
<b>Medallion 16" Gravity Steam Sterilizer - Steam (PMCL P764327-156)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-002
Major Insp. #2 - order PM Pack	P764330-274
<b>Medallion 20" Gravity Steam Sterilizer - Steam (PMCL P764327-156)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 - order PM Pack	P764330-274
<b>Medallion 16" Gravity Steam Sterilizer - Electric w/Chromalox (CHS) (PMCL P764327-156 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-002
Major Insp. #2 for Generator - order PM Pack	P764330-247
Major Insp. #2 for Sterilizer - order PM Pack	P764330-274

\*Descaler 105108 (consumable)



**P.M. PACKS (continued)**

**STEAM STERILIZERS AND STEAM GENERATORS (continued)**

<u>Description</u>	<u>Part Number</u>
<b>Medallion 16" Gravity Steam Sterilizer - Electric w/Electrosteam Solid State (AS) (PMCL P764327-156 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-002
Major Insp. #2 for Generator - order PM Pack	P764330-250
Major Insp. #2 for Sterilizer - order PM Pack	P764330-274
<b>Medallion 16" Gravity Steam Sterilizer - Electric w/Electrosteam Non-Solid State (AS) (PMCL P764327-156 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-246
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-002
Major Insp. #2 for Generator - order PM Pack	P764330-251
Major Insp. #2 for Sterilizer - order PM Pack	P764330-274
<b>Medallion 20" Gravity Steam Sterilizer - Electric w/Chromalox (CHS) (PMCL P764327-156 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 for Generator - order PM Pack	P764330-247
Major Insp. #2 for Sterilizer - order PM Pack	P764330-274

*\*Descaler 105108 (consumable)*

## P.M. PACKS (continued)

### STEAM STERILIZERS AND STEAM GENERATORS (continued)

<u>Description</u>	<u>Part Number</u>
<b>Medallion 20" Gravity Steam Sterilizer - Electric w/Electrosteam Solid State (AS) (PMCL P764327-156 for ster., PMCL P764327-014 for gen.)</b>	
Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-245
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 for Generator - order PM Pack	P764330-250
Major Insp. #2 for Sterilizer - order PM Pack	P764330-274

#### **Medallion 20" Gravity Steam Sterilizer - Electric w/Electrosteam Non-Solid State (AS) (PMCL P764327-156 for ster., PMCL P764327-014 for gen.)**

Major Insp. #1 for Generator - order Descaler* 105108 and PM Pack	P764330-246
Major Insp. #1 for Sterilizer - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 for Generator - order PM Pack	P764330-251
Major Insp. #2 for Sterilizer - order PM Pack	P764330-274

#### **Medallion 24x36" Vacumatic Steam Sterilizer (PMCL P764327-157)**

Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P764324-196
Major Insp. #2 - order PM Pack	P764330-275

#### **Medallion 24x36" Gravity Steam Sterilizer (PMCL P764327-158)**

Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P764324-196
Major Insp. #2 - order PM Pack	P764330-276

## LOW TEMPERATURE STERILIZATION AND EO PRODUCTS

#### **Eagle 3000 20x20" EO Sterilizer (PMCL P764327-007)**

Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 - order PM Pack	P764330-288

\*Descaler 105108 (consumable)

**P.M. PACKS (continued)**

**LOW TEMPERATURE STERILIZATION AND EO PRODUCTS  
(continued)**

<u>Description</u>	<u>Part Number</u>
<b>Eagle 3000 24x36" EO Sterilizer (PMCL P764327-007)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P764324-196
Major Insp. #2 - order PM Pack	P764330-284
<b>Model 5550 EO Analyzer Gas Chromatograph (PMCL P764327-046)</b>	
Major Insp. #1 - order following part(s) separately	
• Cylinder, Permeation	P000039BLN
• Rotor, Cup 10-Port	P000030BLN
• Diaphragm, Sample Pump	P000038BLN
• Filter	NLA
Major Insp. #2 - order following part(s) separately	
• Diaphragm, Sample Pump	P000038BLN
• Filter	NLA
<b>Eagle 2000 20x20" EO Sterilizer (PMCL P764327-067)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 - order PM Pack	P764330-286
<b>Eagle 2000 24x36" EO Sterilizer (PMCL P764327-067)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P764324-196
Major Insp. #2 - order PM Pack	P764330-287
<b>Eagle Series 3017 100% EO Sterilizer - Single/Double Door (PMCL P764327-067)</b>	
Major Insp. #1	PM Pack not required
Major Insp. #2 - order PM Pack	P764330-296
<b>Medallion 16x16" EO Sterilizer (PMCL P764327-076)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-002
Major Insp. #2 - order PM Pack	P764330-283
<b>Medallion 20x20" EO Sterilizer (PMCL P764327-076)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-003
Major Insp. #2 - order PM Pack	P764330-283
<b>Medallion 24x36" EO Sterilizer (PMCL P764327-076)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P764324-196
Major Insp. #2 - order PM Pack	P764330-285
<b>System 1 Processors - Model 89, 90 and 99 (PMCL P764327-251)</b>	
Major Insp. #1 - order PM Pack	200975
Major Insp. #2 - order PM Pack	200975

## P.M. PACKS (continued)

### LOW TEMPERATURE STERILIZATION AND EO PRODUCTS (continued)

<u>Description</u>	<u>Part Number</u>
<b>EO Eagle Monitor - Year 1 (PMCL P764327-100)</b>	
Major Insp. #1 - order following part(s) separately	
• Element, Filter	P129360-802 (4)
• Pump, Sample Carrier	P000019MSI
• Assembly, Charcoal Scrubber	P000069MSI
Minor Insp. #4 - order following part(s) separately	
• Assembly, Charcoal Scrubber	P000069MSI
<b>EO Eagle Monitor - Year 2 (PMCL P764327-100)</b>	
Major Insp. #1 - order following PM Packs	
• PM Pack	P764330-317
• PM Pack	P764330-318
Minor Insp. #4 - order following part(s) separately	
• Assembly, Charcoal Scrubber	P000069MSI
<b>EO Disposer - 50 CFM (PMCL P764327-154)</b>	
Major Insp. #1 - order following part(s) separately	
• Prefilter, Modified	P764326-957
• Kit, Valve Repair/Rebuild	P764326-985
• Kit, Valve Repair/Rebuild	P764327-514
<b>EO Disposer - 100 CFM (PMCL P764327-154)</b>	
Major Insp. #1 - order following part(s) separately	
• Cleaner, Air	P764323-032
• Kit, Solenoid Valve Repair	P764324-031

## DECONTAMINATION

### Reliance 300 Washer

Solenoid/Diverter/Pneumatic Valves	
- order PM Pack	P764328-086

### Reliance 470 Washer

Detergent System - order PM Pack	P764325-795
Solenoid/Diverter/Check Valves	
- order PM Pack	P764325-796
Miscellaneous - order following part(s) separately	
• Spring, Pump	P117950-138
• Pump Cup, Seal Retainer	P117950-139
• Seal, Rubber-Pump	P117950-140
• Seal, Carbon-Pump	P117950-141
• "O" Ring, Pump Cover	P117950-145
• Kit, Repair, 1-3/4" Sol.	P117950-165
• Kit, Repair, 2-1/2" Sol.	P117950-679
• Hose, Neoprene	P117950-690
• Bushing, Spray	P117951-305

**P.M. PACKS (continued)**

**DECONTAMINATION (continued)**

<u>Description</u>	<u>Part Number</u>
<b>Reliance Series 444 Automated Washer (PMCL P764327-257)</b>	
Major Insp. #1 - order PM Pack	P764330-305
Major Insp. #2 - order PM Pack	P764330-306
<b>Reliance Series 120/220 Cart Washer - Single/Double Door (PMCL P764327-009)</b>	
Major Insp. #1 - order PM Pack	P764330-298
Major Insp. #2 - order PM Pack	P764330-299
<b>Reliance Series 430 Washer (PMCL P764327-010)</b>	
Major Insp. #1 - order PM Pack	P764330-303
Major Insp. #2 - order PM Pack	P764330-304
<b>Eagleguard LV Retort Autoclave (PMCL P764327-016)</b>	
Major Insp. #2 - order following part(s) separately	
• Kit, Valve Repair	P764317-687
• Kit, Valve Repair	P764317-688
• Valve, Check	P051770-091
• Valve, Swing Check 1/2" NPT	P150822-354
• Kit, Valve Repair 3/4"	P764078-001
<b>Reliance Series 310 and 980 Washers (PMCL P764327-017)</b>	
Major Insp. #1 - order PM Pack	P764330-300
<b>Reliance Series 410 Washer - Single/Double Door (PMCL P764327-018)</b>	
Major Insp. #1 - order PM Pack	P764330-301
Major Insp. #2 - order PM Pack	P764330-302
<b>Reliance Series 730 Washer (PMCL P764327-019)</b>	
Major Insp. #1 - order PM Pack	P764330-315
Major Insp. #2 - order PM Pack	P764330-316
<b>Reliance Series 777 Automated Washer (PMCL P764327-020)</b>	
Major Insp. #1 - order PM Pack	P764330-307
Major Insp. #2 - order PM Pack	P764330-308
Major Insp. #3 - order PM Pack	P764330-309
<b>Automated Transport System (PMCL P764327-106)</b>	
<i>Replace worn parts as required - frequency depends on usage and specific configuration.</i>	
<b>Sonic Console (PMCL P764327-022)</b>	
Major Insp. #1 - order following part(s) separately	
• Kit, Repair	P752775-091 (3)
• Kit, Valve Repair	P764319-952 (3)
• Kit, Repair	P753448-091 (2)
• Kit, Repair	P756168-091
• Filter, Air	P075261-091

## **P.M. PACKS (continued)**

### **DECONTAMINATION (continued)**

<u>Description</u>	<u>Part Number</u>
<b>Sonic Console w/Special Rinse (PMCL P764327-022)</b>	
Major Insp. #1 - order following part(s) separately	
• Kit, Repair	P752775-091 (3)
• Kit, Valve Repair	P764319-952 (3)
• Kit, Repair	P753448-091 (2)
• Kit, Repair	P756168-091
• Filter, Air	P075261-091
<b>Sonic Bath (PMCL P764327-022)</b>	
Major Insp. #1 - order following part(s) separately	
• Kit, Repair	P752775-091
• Kit, Valve Repair	P764319-952 (2)
• Kit, Repair	P753448-091
• Kit, Repair	P756168-091
<b>3400 Amscomatic Washer/Sterilizer (PMCL P764327-023)</b>	
Major Insp. #1 - order PM Pack	P764330-293
Major Insp. #2 - order PM Pack	P764330-294
<b>Eagle 2000 16x16" PIWS (PMCL P764327-065)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-002
Major Insp. #2 - order PM Pack	P764330-280
<b>Medallion 16x16" PIWS (PMCL P764327-066)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-002
Major Insp. #2 - order PM Pack	P764330-279
<b>Reliance 4800 Automated Washer (PMCL P764327-069)</b>	
Major Insp. #1 - order PM Pack	P764330-311
Major Insp. #2 - order PM Pack	P764330-312
<b>2400 Amscomatic Washer/Sterilizer (PMCL P764327-077)</b>	
Major Insp. #1 - order PM Pack	P764330-291
Major Insp. #2 - order PM Pack	P764330-292
<b>Amscomatic Washer/Sterilizer (PMCL P764327-078)</b>	
Major Insp. #1 - order PM Pack	P764330-289
Major Insp. #2 - order PM Pack	P764330-290
<b>Miele Washer/Disinfector - Model ML200 (PMCL P764327-113)</b>	
Major Insp. #1 - order PM Pack	P000013MLE
<b>Miele Washer/Disinfector - Model ML222 (PMCL P764327-112)</b>	
Major Insp. #1 - order PM Pack	P000012MLE
<b>Eagle 3000 Stage 2 16x16" PIWS (PMCL P764327-153)</b>	
Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-002
Major Insp. #2 - order PM Pack	P764330-281

## **P.M. PACKS (continued)**

### **DECONTAMINATION (continued)**

<u>Description</u>	<u>Part Number</u>
--------------------	--------------------

#### **Eagle 3000 Stage 3 16x16" PIWS (PMCL P764327-153)**

Major Insp. #1 - order Parts Package (for Double Door units - order 2)	P754360-002
---	-------------

Major Insp. #2 - order PM Pack	P764330-282
--------------------------------	-------------

#### **Reliance 333 Washer/Disinfector (PMCL P764327-159)**

Major Insp. #1 - order PM Pack	P764330-149
--------------------------------	-------------

#### **Modulav I and II (PMCL P764327-162)**

Major Insp. #1 - order following part(s) separately

- |                                 |             |
|---------------------------------|-------------|
| • Unit, Hot Quatern             | P759725-001 |
| • Unit, Cold Quatern            | P759725-002 |
| • Kit, Renewal                  | P764316-004 |
| • Sleeve, Rubber                | P764315-998 |
| • Assembly, Diaphragm Operating | P764316-002 |

#### **Reliance 130 Cart and Utensil Washer (PMCL P764327-179)**

Major Insp. #1 - order PM Pack	P764330-265
--------------------------------	-------------

Major Insp. #2 - order PM Pack	P764330-266
--------------------------------	-------------

#### **Reliance Synergy**

Major Insp. #1 - order PM Pack	P764330-818
--------------------------------	-------------

Major Insp. #2 - order PM Pack	P764330-819
--------------------------------	-------------

## **SURGICAL TABLES**

#### **Orthographic II, OrthoVision and OrthoView Orthopedic Fracture Tables (PMCL P764327-042)**

*Replace worn parts as required - frequency depends on usage and specific configuration.*

#### **2080 M, L, RC & IA Surgical Tables (PMCL P764327-043)**

*Replace worn parts as required - frequency depends on usage and specific configuration.*

#### **3080/3085 RC/RL/SP Surgical Tables (PMCL P764327-044)**

*Replace worn parts as required - frequency depends on usage and specific configuration.*

#### **1080 M & L Surgical Tables (PMCL P764327-075)**

*Replace worn parts as required - frequency depends on usage and specific configuration.*

## **SURGICAL LIGHTS**

#### **Quantum Light (PMCL P764327-045)**

*Replace worn parts as required - frequency depends on usage and specific configuration.*

#### **Polaris & Challenge Light (22") (PMCL P764327-047)**

*Replace worn parts as required - frequency depends on usage and specific configuration.*

**P.M. PACKS (continued)**  
**SURGICAL LIGHTS (continued)**

<u>Description</u>	<u>Part Number</u>
<b>Gemini Light (PMCL P764327-049)</b> <i>Replace worn parts as required - frequency depends on usage and specific configuration.</i>	
<b>Examiner 10 Light (PMCL P764327-050)</b> <i>Replace worn parts as required - frequency depends on usage and specific configuration.</i>	
<b>LD202 Light (PMCL P764327-060)</b> <i>Replace worn parts as required - frequency depends on usage and specific configuration.</i>	
<b>SQ-240 Surgical Lighting System (includes SurgiVision) (PMCL P764327-072)</b> <i>Replace worn parts as required - frequency depends on usage and specific configuration.</i>	
<b>Rotoflex and DV22 Series Lights (PMCL P764327-094)</b> <i>Replace worn parts as required - frequency depends on usage and specific configuration.</i>	
<b>SQ-140 Surgical Lighting System (PMCL P764327-117)</b> <i>Replace worn parts as required - frequency depends on usage and specific configuration.</i>	

**SPECIALTY SURGICAL PRODUCTS**

<b>Scrub Sinks (PMCL P764327-051)</b>	
Major Insp. #1 for Unimatic - order PM Pack	P764325-800
Major Insp. #1 for Flexmatic Manual - order following part(s) separately	
• Kit, Water Valve Rebuild	P755715-528
Major Insp. #1 for Flexmatic Infra-red - order following part(s) separately	
• Kit, Water Valve Rebuild	P764319-948
<b>Third &amp; Fourth Generation Steam Warming Cabinets (PMCL P764327-074)</b>	
Major Insp. #1 - order PM Pack	P764330-223
<i>NOTE: Parts not required for Electric Warming Cabinets.</i>	

**PATIENT HANDLING AND TRANSPORT SYSTEMS**

<b>Hausted All Purpose Chairs - Manual and Electric (APC, VIC, ACS, MBC, EPC, ESC) (PMCL P764327-137)</b> <i>Replace worn parts as required - frequency depends on usage and specific configuration.</i>
---



## **P.M. PACKS (continued)**

### **PATIENT HANDLING AND TRANSPORT SYSTEMS (continued)**

<u>Description</u>	<u>Part Number</u>
<b>Hausted Converge, Unicare II &amp; III Stretcher (620, 622, 625, 626, 678, 726, 728, 820, 821, 822, 825, 828, 829 &amp; 878) (PMCL P764327-138)</b> <i>Replace worn parts as required - frequency depends on usage and specific configuration.</i>	
<b>Hausted Horizontal Stretcher Series (416, 428, 432, 436, 466, 471, 472, 492, 493, 496, 4S2, 4S8) (PMCL P764327-139)</b> <i>Replace worn parts as required - frequency depends on usage and specific configuration.</i>	
<b>Hausted Pediatric Stretchers (PD2000, 4160, 6160) (PMCL P764327-140)</b> <i>Replace worn parts as required - frequency depends on usage and specific configuration.</i>	
<b>Hausted Legacy Stretcher Series (PMCL P764327-142)</b> <i>Replace worn parts as required - frequency depends on usage and specific configuration.</i>	
<b>Hausted Gemini Stretcher Series (512, 522, 529 &amp; 562, 578 before 1/31/98) (PMCL P764327-144)</b> <i>Replace worn parts as required - frequency depends on usage and specific configuration.</i>	
<b>Hausted ICU Bed (PMCL P764327-145)</b> <i>Replace worn parts as required - frequency depends on usage and specific configuration.</i>	

### **MISCELLANEOUS EQUIPMENT**

<b>Ultrapure Still (PMCL P764327-087)</b> Major Insp. #1 - order PM Pack	P764330-345
<b>Model C Still (PMCL P764327-087)</b> Major Insp. #1 - order following part(s) separately <ul style="list-style-type: none"><li>• Valve, Swing Check 1/2" NPT</li><li>• Package, Parts</li><li>• Package, Parts - PRV</li></ul>	P150828-442 P754357-006 P754359-002
<b>EcoCycle 10 Processor (PMCL P764327-116)</b> Major Insp. #1 - order PM Pack	P764328-574
<b>SafeCycle 40 (PMCL P764327-121)</b> Major Insp. #1 - order following part(s) separately <ul style="list-style-type: none"><li>• Kit, Valve Repair</li><li>• Kit, Valve Repair</li><li>• Filter, Full Size 1/4" NPT</li><li>• Cartridge, Ametek 1.0 Micron Filter</li></ul> Major Insp. #2 - order following part(s) separately <ul style="list-style-type: none"><li>• Cartridge, Ametek 1.0 Micron Filter</li></ul>	P764328-644 (6) P764328-638 (3) 400727 (2) 400688 400688

## LAMPS

LAMP - FOR 18" LONG BJRS, BJRA, BM AND AJR SERIES.  
DOOR GLASS USED TO SEAL BOTTOM OF REFLECTOR

**NLA** LAMP, 150W, 230V, 150P25/10SB, MED. BASE

LAMP - FOR 22" ROTOFLEX SERIES - C-22

**P045704-091** LAMP, 150W, 120V, 150P25/2SB, MED. BASE

LAMPS - FOR 22" SERIES SURGICAL LIGHTS - C-22-O, C-22-S, CD-22-O, CD-22-S, C-22-E, C-22-T AND DV22 SERIES. BOTTOM OF REFLECTOR OPEN - NO DOOR GLASS USED (NOTE: LAMP FOR PX22 - NLA)

**NLA** LAMP, 150W, 24V, 150P25/10, MED. BASE

**NLA** LAMP, 150W, 230V, 150P25/10SB, MED. BASE

**NLA** LAMP, 150W, 120V, 150P25/8SB, MED. BASE

LAMPS FOR CHALLENGE 22" AND POLARIS LIGHT

**P129249-001** LAMP, 235W, 33V.

**P764319-682** LAMP, BOXED (12 EACH P-129249-001 PER BOX)

LAMPS FOR EXAMINER 10 LIGHT

**P084079-033** LAMP, 120V, 150W (DOMESTIC)

**P129352-157** LAMP, 33V, 235W (INTERNATIONAL)

LAMP FOR MODEL 180 LIGHT

**P129357-106** LAMP, 33V, 235W

LAMP FOR POLARIS SL 34" AND MODEL 900 EXAM LIGHT

**P093902-608** LAMP, 22V, 50W

LAMP FOR GEMINI 22" LIGHT

**P129359-934** LAMP, 22V, 200W

LAMP FOR DMC 34" LIGHT

**P764324-049** LAMP, 2-WAY, 115V, 150/400W

LAMP FOR QUANTUM LIGHT, LD202 AND SQ240

**P129362-228** LAMP, 22V, 220W

LAMP FOR HARMONY

**P093926-047** LAMP, 20V, 180W

LAMP FOR DEEPSITE

**P136818-113** LAMP, XENON ASS'Y

## **STERILIZABLE/DISPOSABLE SURGICAL LIGHT HANDLES**

- P027183-033** HANDLE, STERILIZABLE (ROTOFLEX)
- P093184-001** HANDLE, STERILIZABLE (FOR CHALLENGE 22, POLARIS, SQ240)
- P146655-001** HANDLE, STERILIZABLE (QUANTUM)
- P136820-001** HANDLE, STANDARD (HARMONY)
- P093926-125 COVER, LIGHTHANDLE (STERILE, DUAL PACK)
  - P093926-328 COVER, LIGHTHANDLE (STERILE, SINGLE PACK)
- P136818-103** SLEEVE, DISPOSABLE STERILE (DEEPSITE)

## 2000 SERIES STERILIZERS W/PRINTCON P/C BOARDS

Old Part Number*	Where Used	New Part Number
P146654-002 (NLA)	2300 w/Olivetti Printer	P762650-481
P146654-003 (NLA)	Gas Sterilizer w/Olivetti Printer and Air Flow Alarm Upgrade	P762650-482
P146654-009 (NLA)	Steam Sterilizers, PIWS and pre-1987 Gas without Air Flow Alarm Upgrade w/Olivetti Printer	P762650-483
P146654-012 (NLA)	2400 w/Olivetti Printer	P762650-484
P146654-016 (NLA)	Steam Sterilizers, PIWS and pre-1987 Gas without Air Flow Alarm Upgrade w/Epson Printer	P762650-485
P146654-017 (NLA)	Gas Sterilizer w/Epson Printer and Air Flow Alarm Upgrade	P762650-486
P146654-021 (NLA)	2300 w/Epson Printer	P762650-487
P146654-022 (NLA)	2300 Floorloader w/Epson Printer	P762650-488
P146654-023 (NLA)	2400 w/Epson Printer	P762650-489
P146654-024 (NLA)	Terminal Sterilizer w/Olivetti Printer	P762650-498

**\* Old Part Numbers are No Longer Available; sub to New Part Numbers.**

## PRINTED CIRCUIT BOARDS

### Eagle 2000/2400 Pulstar 2000/3000/5000

<b>New</b>		
<b>Part Number</b>		<b>Description</b>
P141198-310 .....		Power Supply, 5 Fuse Printcon
P146645-226 .....	See Note 1	Panel, Primary Control
P146302-001 .....		Panel, Primary Control W/S
P764320-963 .....	See Note 2	Power Supply I.R.C.
P146517-001 .....		Panel, N.O.E. Primary Control
P146528-001 .....		Panel, N.O.E. Primary Control PIWS
P146534-001 .....		Panel, N.O.E. Primary Control VAC
P146586-001 .....		PC Board #1
P146586-002 .....		PC Board #1 Pulstar 2000/3000 (Before 2/84)
P146586-003 .....		PC Board #1 Eagle PIWS
P146586-004 .....		PC Board #1 Pulstar 5000
P146586-006 .....		PC Board #1 Eagle
NLA .....		PC Board #1
P146586-008 .....		PC Board #1 Pulstar 2000/3000 (After 2/84)
P146586-101 .....	See Note 3	PC Board #1 Eagle 2400
P146586-101 .....	See Note 3	PC Board #1 Eagle 2400
P146586-101 .....	See Note 3	PC Board #1 Eagle 2400
P146586-082 .....		PC Board #1 RV ETO Gas
P146586-101 .....		PC Board #1 AWCC 2400
NLA .....		PC Board #2 Eagle
P146588-002 .....		PC Board #2 Pulstar 2000/3000 (Before 2/84)
P146588-004 .....		PC Board #2 Pulstar 5000
NLA .....		PC Board #2 Eagle
P146588-007 .....		PC Board #2 Eagle
P146588-008 .....		PC Board #2 Pulstar 2000/3000 (After 2/84)
P146588-101 .....	See Note 3	PC Board #2 Eagle 2400
NLA .....		PC Board #2 Eagle
P146588-058 .....		PC Board #2 Eagle
P146588-101 .....	See Note 3	PC Board #2 Eagle
P146588-101 .....	See Note 3	PC Board #2 Eagle 2400
P146588-101 .....		PC Board #2 AWCC 2400
P146590-001 .....		PC Board #3 Eagle, Pulstar
P146590-007 .....		PC Board #3 Eagle
P146590-054 .....		PC Board #3 Eagle
P146590-058 .....		PC Board #3 Eagle
P146590-101 .....	See Note 3	PC Board #3 2400
P146590-101 .....		PC Board #3 AWCC 2400
P146633-031 .....		Motherboard 2200 & 2400
P146633-081 .....		PC Board Printcon Display
P141198-310 .....	See Note 4	Power Supply Printcon
P146645-226 .....		Panel Primary Control
P146645-227 .....		Panel, N.O.E. Primary Control Gas/Aer.
P146648-001 .....		PC Board #4 Eagle 2400
P762650-483 .....	See Note 5	PC Board Printcon 2000
P762650-484 .....	See Note 5	PC Board Printcon 2400

P150475-878 .....	Panel Primary
NLA .....	Motherboard with Chassis, Gravity Units
NLA .....	Motherboard with Chassis, All Except 2200 & 2400
NLA .....	Motherboard and Chassis PIWS
P134471-039 .....	Panel Primary Control

#### **Eagle 2300**

P146586-051 .....	PC Board #1 Eagle 2300
P146588-051 .....	PC Board #2 Eagle 2300

#### **Pulstar 1200**

NLA .....	PC Board #1
NLA .....	PC Board #2

#### **180 Surgical Lights**

P146645-639 .....	PC Board 180 Light
-------------------	--------------------

#### **180 Sterilizers**

P093904-600 .....	PC Board 180 Control
P146645-682 .....	PC Board Printer/Display

#### **900 Sterilizer**

P093904-700 .....	PC Board Revision A
P093904-702 .....	PC Board Revision B

#### **Eagle 3000 (Steam Sterilizer) Stage II**

P146653-204 .....	PC Board Printer
P146653-201 .....	PC Board Control (See Note 6)

#### **Eagle 3000 (Gas Sterilizer) Stage III**

P146655-356 .....	Control Board
P093910-531 .....	Printer Assembly

#### **Eagle 3000 (Vac Floor Loading Sterilizer)**

P146653-204 .....	PC Board Printer/Display
P146653-201 .....	PC Board Control

#### **Eagle 2300 (Vac Floor Loading Sterilizer)**

P146586-070 .....	PC Board #1
P146588-070 .....	PC Board #2
P146590-070 .....	PC Board #3

#### **Eagle 3400**

P146655-356 .....	PC Board Control
P146655-598 .....	See Note 7 .... PC Board Interface
P093910-531 .....	Printer Assembly
P093910-647 .....	See Note 7 .... Eprom

#### **AWCC 2400**

P146586-101 .....	PC Board #1
P146588-101 .....	PC Board #2
P146590-101 .....	PC Board #3
P146648-001 .....	PC Board #4

**Century/Millennium 116/120**  
**using P146656-185 Control Board**

P146659-005 .....	I/O Board
P146659-070 .....	Connector Board
P146656-182 .....	Printer Board

**Century/Millennium 116/120**  
**using P146659-065 & P146659-069 Control Boards**

P146659-066 .....	I/O Board (EN285)
P146659-070 .....	Connector Board
P146656-182 .....	Printer Board

**AMSCO Century Medium Sterilizer (26 x 37.5 Chamber)**

P141214-630 .....	PC Board, Interface
P146656-182 .....	PC Board, Printer
P141215-202 .....	PC Board, CPU
P146659-008 .....	PC Board, Input/Output #1
P146659-005 .....	PC Board, Input/Output #2
P136810-006 .....	Power Supply

---

**Note 1:** P146163-001(NLA) replaced by P146645-226 Panel Primary Control.

**Note 2:** P146402-001(NLA) replaced by P764320-963 Kit Power Supply Conversion.

**Note 3:** The following P.C. Boards are no longer available. Replacement boards are as indicated:

**Eagle 2400 Washer Sterilizer**

P146586-067(NLA) replaced by P146586-101 PC Board #1  
P146588-067(NLA) replaced by P146588-101 PC Board #2  
P146586-048(NLA) replaced by P146586-101 PC Board #1  
P146588-048(NLA) replaced by P146588-101 PC Board #2

**AWCC 2400**

P146586-078(NLA) replaced by P146586-101 PC Board #1  
P146588-078(NLA) replaced by P146588-101 PC Board #2  
P146590-078(NLA) replaced by P146590-101 PC Board #3  
(continue to use P146648-001 as PC Board #4)

**Note 4:** Power supply P146645-042 (NLA) is replaced by P141198-310. This is a direct replacement for all applications except Eagle 2014 Washer Sterilizer. The Eagle 2014 Washer Sterilizer requires power supply and relay interface kit P150822-969.

**Note 5:** P146651-001(NLA) replaced by P762650-483 PC Board Printcon 2000. P146651-032(NLA) replaced by P762650-484 PC Board Printcon 2400.

**Note 6:** P146653-201 Control Board does not include Eprom chip. Refer to Eprom for part number.

**Note 7:** P146655-366(NLA) Interface P.C. Board Assembly replaced by P146655-598. P056400-091 replaced by P093910-647 programmed Eprom.

---

**Important Note**

**Per STERIS's new return goods policy, printed circuit boards are not eligible for return and restocking.**

## EAGLE 3000 STAGE 3 CONTROLS

<u>Description</u>	<u>Part Number</u>
<b>Control: Standard Vertical Column</b>	
Power Supply (2 board control) .....	P136810-006
Power Supply (single board control) .....	P093910-541
PC (Printed Circuit) Board, Control (w/2K RAM) .....	P146655-356
• EPROM* (Rev. X) .....	P093913-XXX
• RAM (2K) .....	P150828-625
• Triac .....	P129360-770
PC (Printed Circuit) Board, Control (w/8K RAM**) .....	P146655-819
• RAM (8K) .....	P129362-145
PC (Printed Circuit) Board Connector .....	P146655-357
• Fuse, 3 amp .....	P764317-778
Touch Pad Panel .....	P136806-685
Motor, Paper Take-up .....	P093909-399
Fuse, Line Filter (5 amp) .....	P764317-463
Line Filter .....	P093910-584
<b>If Thermal Printer Equipped</b>	
Printer Assembly .....	P093910-531
• Printer Only .....	P755716-005
Paper (5 Rolls) .....	P129359-008
Control Interface Board .....	P146655-598
<b>If Impact Printer Equipped</b>	
<b>Printer Assembly Before 1-96</b>	
Printer Board Assembly, Impact .....	P338517-170
• Printer Assembly w/Cable .....	P338517-173
• Board Only (Epson Type) .....	P338517-164
• Cartridge, Ink Printer .....	P150828-440
• Paper, Single Ply, 3 Rolls .....	P764324-057
• Paper, Double Ply, 3 Rolls .....	P150828-420
<b>Printer Assembly After 1-96</b>	
Printer Board Assembly, Impact .....	P093911-498
• Printer Only .....	P093918-051
• Board Only .....	P146656-182
• Cartridge, Ink Ribbon .....	P150828-440
• Paper, Single Ply, 3 Rolls .....	P764324-057
• Paper, Double Ply, 3 Rolls .....	P150828-420
<b>3000 SL Impact Printer</b>	
Printer Board Assembly .....	P093911-459
• Printer Only .....	P093911-502
• Board Only .....	P056401-989
• Cartridge, Ink Ribbon .....	P150828-440
• Paper, Single Ply, 3 Rolls .....	P129362-819

*\*NOTE: When ordering P093913-XXX, Rev. X, EPROMs refer to the part number on the label in the control column. If one control board is present, order one EPROM. If two control boards are present, order two EPROMs (part numbers may be different so check both EPROMs on the control boards).*

*\*\*NOTE: The 8K RAM will fill all the pins in the socket (28 pins). The 2K RAM has only 24 pins.*



## PC Boards Return List

<b>Description</b>	<b>Part Number</b>
<b>System 1 Model 89/90 Units:</b>	
CPU Board .....	300086
Printer Board .....	300106
<b>System 1 Model 99 Units:</b>	
Interface Board .....	P093915-200
CPU Board .....	P146665-002
<b>Century Sterilizers (Small and Medium):</b>	
CPU Board .....	P146659-065
CPU Board .....	P146659-069
CPU Board .....	P141215-202
I/O Board .....	P146659-005
I/O Board .....	P146659-008
I/O Board .....	P146659-066
Display Interface Board .....	P141214-630
Printer Board .....	NLA
<b>Eagle 3000 Stage 3 Sterilizers and Washers:</b>	
CPU Board .....	P146655-356
CPU Board .....	P146655-819
<b>333 Washer:</b>	
CPU Board .....	P117910-844
I/O Board .....	P117909-442
<b>Harmony Lighting System:</b>	
Wall Unit Controller .....	P146667-035
Lamp Head Controller Board .....	P146667-106
Standard Power Status Board .....	P136820-101
HUB Interface Board .....	P136820-105
Analog Video Interface Board .....	P136820-107
DeepSite Interface Board .....	P136820-109
Backup Power Status Board .....	P136820-103
Small Lamp Head Controller Board - LA 300 .....	P136820-111
<b>Eagle 2000 and Eagle 3000:</b>	
P/C Assembly Basic .....	P146586-001
P/C Assembly Basic Instrument Wash .....	P146586-003
P/C Board .....	NLA
Programmed Board #1 .....	P146586-101
P/C Card Assembly, Vac and Gas .....	NLA
P/C Card Assembly, Vac and Gas .....	NLA
P/C Board Assembly .....	P146588-007
P/C Card Assembly .....	NLA
P/C Board #3 Eagle .....	P146590-058
Control Board - Eagle 3000 .....	P146655-356
Control Board - Eagle 3000 .....	P146655-819
Supply-Power For Sgl Board .....	P093910-541
Supply-Power For Dbl Board .....	P136810-006
Control Board: Eagle 3000 .....	P146653-201
Printer Board .....	P146653-204

## TABLE LIFT CYLINDER

Table	New Lift Cylinders
1080 .....	P136406-001
2080 .....	P136375-001

## COMMON ACCESSORIES

Clark Socket .....	P043348-091
Knob/Socket .....	P046145-091
Armboard only .....	P056130-001
• Bracket Repair Kit .....	P764316-574
Large Post Kidney .....	P626397-001
Small Post Kidney .....	P626397-002
Traction Boot .....	P136804-149

## ORTHOGRAPHIC 2

T-Shaped Sacral Rest Assembly* .....	P136804-138
Wedge Shaped Sacral Rest Assembly* .....	P136804-139
Child Sacral Rest Assembly* .....	P056397-105
Post Assembly, 11" * .....	P129186-109
Post Assembly, 9" * .....	P755716-035
Child's Post Assembly * .....	P056397-117

## ORTHOVISION

T-shaped Sacral Rest .....	P134469-365
• Pad .....	P134469-364
Sacral Rest Assembly .....	P134469-225
• Pad .....	P150830-170
Post Assembly .....	P093910-572

*\* All three posts and sacral rests are interchangeable.*

### 3080/3085 HYDRAULIC CYLINDERS

**3080 "04" Hydraulic Cylinders** - All S/N starting with "04", the original hydraulic cylinders are no longer available. Retrofit kits are listed below with the cylinder p/n's for replacement after units have been updated with retrofit kits.

*Note: p/c - one kit needed per cylinder; p/t - one kit per table, these kits include right & left cylinders, if applicable.*

Articulation	Original "04" P/N (qty)	Kit to Retrofit	Cylinder after Retrofit
Floor Lock Cylinders	P146653-540 (2); P146653-777 (1) - NLA	P150823-963 p/c	P056397-275 (3)
Seat Cylinder	P141210-098(1) - NLA	P150823-983 p/t	P134469-139(1)
Side Tilt Cylinder	P141210-107 (1)	N/A	N/A
Back Section Cylinders	P141210-096 (1); P141210-099 (1) - NLA	P150823-982 p/t	P134469-131 (2)
Leg Section Cylinder Left	P146655-463 (1) - NLA	P150823-981 p/t	P141210-273 L.H. (1)
Leg Section Cylinder Right	P146655-464 (1) - NLA	P150823-981 p/t	P141210-274 R.H. (1)

**3080/3085 Bieri "B4" Hydraulic Cylinders** - There are 3 different hydraulic configurations on the 3080/3085 Tables that the S/N begins with "B4". Note the S/N of the unit before determining cylinder replacement p/n.

Articulation	Bieri 1 - S/N up to B426594-084 and S/N B432195-001 thru B432195-100	Bieri 2 - S/N from B426594-085 up to B403797-027 (with exceptions noted in Bieri 1)	Bieri 3 - S/N from B403797-028 to present; and all 3085 Tables
Floor Lock Cylinder #1	P056397-274	P056397-274	P056397-274
Floor Lock Cylinder #2	P056397-271	P056397-275	P056397-275
Floor Lock Cylinder #3	P056397-272	P056397-275	P056397-275
Seat Cylinder	P134469-142	P134469-133	P134469-133
Side Tilt Cylinder	P141210-108	P141210-108	P141210-108
Back Section Cylinder Left	P134469-140	P134469-131	P134469-131
Back Section Cylinder Right	P134469-145 (Kit to install P134469-132)	P134469-132	P134469-132
Leg Section Cylinder Left	P141210-145	P141210-145	P141210-145
Leg Section Cylinder Right	P141210-146	P141210-146	P141210-146

## HAUSTED CYLINDER ASSEMBLIES, COLUMNS & MEDIJACKS

<u>Description</u>	<u>Part Number</u>
<b>Hausted Hydraulic Cylinder:</b>	
Unicare Series .....	065087
Gemini .....	067674
Gemini 578 bases (pre 95) .....	065087
Converge (pre 01/95) .....	-Hd. End-061517
Converge (pre 01/95) .....	-Ft. End-061520
2700 Bed .....	-Hd. End-061517
2700 Bed .....	-Ft. End-061520
<b>Linear Column:</b>	
EPS/ESC Chairs .....	068972
Converge (post 01/95) .....	066516
Vision 2000 (793 Base) .....	066212
Legacy (893 Base) .....	066129
Legacy (892 Base) .....	-Hd. End-066109
Legacy (892 Base) .....	-Ft. End-066107
Horizon Electric (4E2, 5E8) .....	068983
<b>Medijack Self-contained:</b>	
APC/VIC/MBC Chairs .....	061122
Horizon (432, 462, 492, 578 after 01/95) .....	-Hd. End-068209
Horizon (432, 462, 492, 578 after 01/95) .....	-Ft. End-068207
Horizon (4160, 468, 471, 472, 578 after 02/98) .....	-Hd. End-075309
Horizon (4160, 468, 471, 472, 578 after 02/98) .....	-Ft. End-075310

## SURGICAL AND ACCESSORY PADS

Notes:

1. Only **current** capital equipment *Table Top Pads* are available in the 2" Foam and the 2-1/4" TLT. All other *Table Top Pads* and *Accessory Pads* are **only** available in the TLT offering.
2. *Foam* - 2" Thick basic foam pads.
3. *TLT* (Tri Layer Technology) - are our standard 2-1/4" pressure management offering.
4. *N/A* - this option is not offered for that style pad.
5. *NLA* - this pad is no longer offered and/or available.

### 1063/1070/1080/2080 SURGICAL TABLE PADS

<u>Pad Description:</u>	<u>Foam P/N:</u>	<u>TLT, 2-1/4" P/N:</u> (Unless otherwise noted)
Pad, Complete Foam Set, 2"	NLA	N/A
• Pad, Leg, 2"	P093074-001	N/A
• Pad, Back/Seat, 2"	P093075-001	N/A
• Pad, Head, 2"	P093076-001	N/A
Pad, Complete "TLT" Set, 2-1/4"	N/A	P150830-177
• Pad, Head, 2080 "TLT" 2-1/4"	N/A	P150830-166
• Back/Seat, 2080 "TLT" 2-1/4"	N/A	P150830-165
• Pad, Leg, 2080 "TLT" 2-1/4"	N/A	P150830-164
Pad, Foot Section for Hess Drain Tray, 2-1/4"	N/A	P755716-780
Pad, Foot Extension, 1"	N/A	P150830-185
Pad, Foot Extension, 2"	N/A	P150830-185
Pad, Head Section, 3-1/4"	N/A	P617101-164
Pad, Back Section, 3-1/4"	N/A	P617101-165
Pad, Leg Section, 3-1/4"	N/A	P617101-166
Pad, Head Section, 4-1/4"	N/A	P616978-001
Pad, Back Section, 4-1/4"	N/A	P616979-001
Pad, Leg Section, 4-1/4"	N/A	P616980-001
Pad, Mattress, DAD, 2-1/4" x 73", (Covers 2080 Tabletop)	N/A	NLA
Pad, Mattress, DAD, 2-1/4" x 85" (Covers 2080 Tabletop with Foot Section)	N/A	P150830-165

### OLD STYLE 1080 TABLE PADS

Pad, Complete Foam Set, 2-1/4"	N/A	P755716-781
• Pad, Head Section, 2-1/4"	N/A	P755716-782
• Pad, Back Section, 2-1/4"	N/A	P755716-783
• Pad, Foot Section, 2-1/4"	N/A	P755716-784

### 3080/3085 SURGICAL TABLE PADS

Pad, Set, 2-1/4"	N/A	P150830-176
• Pad, Head, 2-1/4"	N/A	P150830-175
• Pad, Back/Seat, 2-1/4"	N/A	P150830-173
• Pad, Leg Section, 2-1/4"	N/A	P150830-174
Pad, Complete, 2"	P129360-589	N/A
• Pad, Leg Section, 2"	P136806-811	N/A
• Pad, Head Section, 2"	P136806-813	N/A
• Pad, Back/Seat Section, 2"	P136806-815	N/A
Pad, Image Intensification Top 1" thick x 86" long	N/A	P150830-194
Pad, Image Intensification Top 2" thick x 62" long	N/A	P150830-211
Pad, Lateral Brace Roll	N/A	P093910-811
Pad, 4-1/4", Leg Transfer Board	N/A	P150830-189
Pad, Lightweight Patient Transfer Board	N/A	P150830-188
Pad, Infant Surgery Attachment	N/A	P755716-785
Pad, Hourglass Arm & Hand Table	N/A	P134469-229
Pad, Neurosurgeon Armrest	N/A	P755716-811

### 3080/3085 SURGICAL TABLE PADS (continued)

<u>Pad Description:</u>	<u>Foam P/N:</u>	<u>TLT, 2-1/4" P/N:</u> (Unless otherwise noted)
Pad, Lightweight Patient Transfer Board	N/A	P150830-187
Pad, Infant Table	N/A	P136807-025
Pad, Image Intensification, Top (old style) 2-1/4" thick x 74" long	N/A	P755716-786
Pad, I.A. Extender, Tabletop Ext.	N/A	P150830-193
Pad, Fem/Pop Board	N/A	P150830-186
Pad, 3080 Leg Section Wide X-ray Top	N/A	P134469-097
Pad, 3080 Head Section Wide X-ray Top	N/A	P134469-098
Pad, 3080 Back/seat Section, Wide X-ray Top	N/A	P134469-099
Pad, Ophthalmic Headrest	N/A	P056124-001
Pad, Cone Shaped Headrest	N/A	P150830-195
Pad, Cone Shaped Headrest	N/A	P150830-196
Pad, Captains Chair Main Body	N/A	P056397-547
Pad, Captains Chair Right-hand Shoulder	N/A	P056397-548
Pad, Captains Chair Left-hand Shoulder	N/A	P056397-548
Pad, Captains Chair Headrest	N/A	P056397-540
Pad, Morgan Pelvic Foot Extension	N/A	P080441-091
Pad, Morgan Pelvic Rest, 4-1/4" thick - Back/Seat Section w/larger perineal cutout	N/A	P080524-091
Pad, Morgan Pelvic Rest, 4-1/4" thick	N/A	P097021-091
Pad, Foot Extension, 1" thick	N/A	P150830-185
Pad, Bierhoff Legholder, 1"	N/A	P755716-787
Pad, Nager Headrest	N/A	P755716-788
Pad, 3080 I.A. Extender	N/A	P150830-193
Pad, Foot Extension	N/A	P150830-185
Pad, Culdoscopy Knee Rest	N/A	P150830-183
Pad, 3080 Neuro Back Extension	N/A	P136807-094
Pad, Neuro Elbow Rest	N/A	P755716-811
Pad, 2080 Perineal Filler	N/A	P617101-217
Pad, 2080 Perineal Filler, 3"	N/A	P755716-065
Pad, 3080 Perineal Filler	N/A	P617101-698
Pad, Disposable Foam Insert for Multipurpose Legholders	NO-002	N/A

### ORTHOGRAPHIC 2 AND ORTHOVISION ORTHOPEDIC TABLE PADS

Pad, Table Top, Foam, 2"	P136804-103	N/A
Pad, Table Top, DAD, 2-1/4"	N/A	P150830-172

### ORTHOPEDIC ACCESSORY PADS

Pad, Ortho II/Vision/View Set	N/A	P150830-234
• Pad, Table Top	N/A	P150830-172
• Pad, T-shape Sacral Rest	N/A	P150830-169
• Pad, Wedge Shape	N/A	P150830-170
• Pad, I.M. Nailing Sacral Rest (Hip Rest) - 2-1/4"	N/A	P150830-171
Pad, Patient Transfer Board, Old Style	N/A	P755716-789
Pad, Patient Transfer Board, Old Style	N/A	P755716-793
Pad, Lightweight Patient Transfer Board	N/A	P150830-187
Pad, Image Amplification Board	N/A	P150830-191
Pad, Image Amplification Board, Old Style	N/A	P150830-191
Pad, Arm & Hand Table, 2-1/4" (BF00-026)	N/A	P755716-790
Pad, Arm & Hand Table, 3" (BF00-026)	N/A	P755716-791
Pad, Arm & Hand Table, 4" (BF00-026)	N/A	P755716-792
Boot, Traction	P136804-149	N/A
Pad, Universal Legholder, Molded	P136804-040	N/A
Pad, Roll, Popliteal Support	P093898-119	N/A
Pad, Disposable Perineal Support	NO-001	N/A
Pad, Leg Rest	N/A	BF00062
Pad, Hip Spica	N/A	P056397-121
Pad, Child Sacral Rest, Foam	N/A	P150830-190

## ORTHOPEDIC ACCESSORY PADS (continued)

<u>Pad Description:</u>	<u>Foam P/N:</u>	<u>TLT, 2-1/4" P/N:</u> (Unless otherwise noted)
Pad, Tibia Countertraction	N/A	P136807-115
Pad, Arm Extension Assembly	N/A	NLA
Pad, Knee Support	N/A	NLA
Pad, 2" Board	N/A	P150830-191

### ARMBOARD PADS

Pad, 2-1/4"	N/A	P150830-167
Pad, 3-1/4" thick	N/A	P150830-168
Pad, 4-1/4" thick	N/A	P617101-098
Pad, Multi Task	N/A	P150830-178
Pad, 1" (for Surgigraphic only)	P200050-528	N/A

### DELIVERY 4 TABLE PADS

Pad, Seat/back Section, 2-1/4"	N/A	NLA
Pad, Legboard (16x29x2-1/4")	N/A	P755716-796
Pad, Legboard (24x29x2-1/4")	N/A	P755716-797

### ORTHOGRAPHIC (FUCHES) TABLE PADS

Pad, Lateral Support Assembly (P142781-111)	NLA	NLA
Pad, Dorsal Plate Assembly (NLA)	NLA	NLA
Pad, Head Rest-Seat and Head Sections w/Irons (P142773-100)	NLA	NLA
Pad, Seat Section - Seat and Head Sections w/Irons (P142773-101)	NLA	NLA
Pad, Leg Section - Leg Plate Assembly (NLA)	NLA	NLA
Pad, Pelvis Plate - Nailing Device (P142781-103)	NLA	NLA
Pad, Roll - Nailing Device (NLA)	NLA	NLA

### ARNOLD ENDO-URO 1100 TABLE PADS

Pad, Back Section, 2-1/4"	N/A	P755716-798
Pad, Legboard, 2-1/4" (NLA)	NLA	NLA
Pad, Legboard, 2-1/4" x 35" long	N/A	P755716-799
Pad, Legboard, 2-1/4" x 25" long	N/A	P755716-800

### ARNOLD ENDOGRAPHIC TABLE PADS

Pad, Back Section (NLA)	NLA	NLA*
Pad, Leg Section (NLA)	NLA	NLA*

### GU110 TABLE PADS

Pad, Back Section, 2-1/4" (P045977-091)	NLA	NLA
Pad, Foot Section, 2-1/4" (NLA)	NLA	NLA
Pad, Back Section, 3" (NLA)	NLA	NLA
Pad, Foot Section, 3" (NLA)	NLA	NLA

### 500 TABLE PADS

Pad, Head Section, 2-1/4"	NLA	P755716-805
Pad, Foot Section, 2-1/4"	NLA	P755716-806
Pad, Fitzhugh Backrest (NLA)	NLA	NLA

\* Reupolster only

### 800 TABLE PADS

<u>Pad Description:</u>	<u>Foam P/N:</u>	<u>TLT, 2-1/4" P/N:</u> (Unless otherwise noted)
Pad, Back Section, 2-1/4" (NLA)	NLA	NLA
Pad, Foot Section, 2-1/4" (NLA)	NLA	NLA

### 800W TABLE PADS

Pad, Back Section, 2-1/4" (NLA)	NLA	NLA
Pad, Foot Section, 2-1/4" (NLA)	NLA	NLA



## BRISTOL TO U.E. RECORDER/CONTROLLER REPLACEMENT CHART

### Eagle Small & Medium Gravity/Isothermal and Vacamatic Sterilizers:

Existing Bristol Recorder/Controller .....	P764321-964
U.E. Recorder Replacement Kit .....	P764321-964
U.E. Recorder - Available Spare Parts:	
· U.E. Charts (Box of 100) °F .....	P056396-503
· U.E. Charts (Box of 100) °C .....	P056396-506
· U.E. Pen Arm Asm .....	P764321-968
· U.E. Pens .....	P764321-741
· U.E. Thermal Asm .....	P764321-969
· U.E. TS-1 Switch TS-2 .....	P764321-972
· U.E. Chart Drive (60 Cycle) .....	P764321-967
· U.E. Chart Drive (50 Cycle) .....	P764322-096
· U.E. Recorder/Controller .....	P141198-270

### Eagle Small Washer/Sterilizers:

Existing Bristol Recorder/Controller .....	P764321-965
U.E. Recorder Replacement Kit .....	P764321-965
U.E. Recorder - Available Spare Parts:	
· U.E. Charts (Box of 100) °F .....	P056396-763
· U.E. Charts (Box of 100) °C .....	P056396-764
· U.E. Pen Arm Asm .....	P764321-968
· U.E. Pens .....	P764321-741
· U.E. Thermal Asm .....	P764321-970
· U.E. TS-1 Switch TS-2 .....	P764321-972
· U.E. Chart Drive (60 Cycle) .....	P764321-967
· U.E. Chart Drive (50 Cycle) .....	P764322-096
· U.E. Recorder/Controller .....	P141198-271

### Medallion Small & Medium G.P./Lab Sterilizers - Medallion Medium Vac A/S Sterilizers Manufactured After 1971 - Medallion Medium Vac S Sterilizers Manufactured Before 1971 - All Medallion X Sterilizers:

Existing Bristol Recorder/Controller	
(2 switch) .....	P764321-997
(3 switch) .....	P764321-998
U.E. Recorder Replacement Kit .....	P764321-997
U.E. Recorder - Available Spare Parts:	
· U.E. Charts (Box of 100) °F .....	P056396-503
· U.E. Charts (Box of 100) °C .....	P056396-506
· U.E. Pen Arm Asm .....	P764321-968
· U.E. Pens .....	P764321-741
· U.E. Thermal Asm .....	P764321-969
· U.E. TS-1 Switch TS-2 .....	P764321-972
· U.E. Chart Drive (60 Cycle) .....	P764321-967
· U.E. Chart Drive (50 Cycle) .....	P764322-096
· U.E. Recorder/Controller .....	P760265-561

### Medallion Small & Medium Lab/Iso Sterilizers:

Existing Bristol Recorder/Controller	
(2 switch) .....	P764321-997
(3 switch) .....	P764321-998
U.E. Recorder Replacement Kit .....	P764321-998
U.E. Recorder - Available Spare Parts:	
· U.E. Charts (Box of 100) °F .....	P056396-503
· U.E. Charts (Box of 100) °C .....	P056396-506
· U.E. Pen Arm Asm .....	P764321-968
· U.E. Pens .....	P764321-741
· U.E. Thermal Asm .....	P764321-969
· U.E. TS-1 Switch TS-2 .....	P764321-972
· U.E. Chart Drive (60 Cycle) .....	P764321-967
· U.E. Chart Drive (50 Cycle) .....	P764322-096
· U.E. Recorder/Controller .....	P760265-562

# BRISTOL TO U.E. RECORDER/CONTROLLER REPLACEMENT CHART (Cont'd.)

## Medallion Small Washer/Sterilizers:

Existing Bristol Recorder/Controller .....	P764321-999
U.E. Recorder Replacement Kit .....	P764321-999
U.E. Recorder - Available Spare Parts:	
· U.E. Charts (Box of 100) °F .....	P056396-763
· U.E. Charts (Box of 100) °C .....	P056396-764
· U.E. Pen Arm Asm .....	P764321-968
· U.E. Pens .....	P764321-741
· U.E. Thermal Asm .....	P764321-970
· U.E. TS-1 Switch TS-2 .....	P764321-972
· U.E. Chart Drive (60 Cycle) .....	P764321-967
· U.E. Chart Drive (50 Cycle) .....	P764322-096
· U.E. Recorder/Controller .....	P760265-565

## Medallion Small & Medium EO (Gas) Sterilizers:

Existing Bristol Recorder/Controller .....	P764322-002
U.E. Recorder Replacement Kit .....	P764322-002
U.E. Recorder - Available Spare Parts:	
· U.E. Charts (Box of 100) °F .....	P755715-857
· U.E. Charts (Box of 100) °C .....	N/A
· U.E. Pen Arm Asm .....	N/A
· U.E. Pens (Red) .....	P764321-958
· U.E. Pens (Blue) .....	P764321-742
· U.E. Thermal Asm .....	P764321-971
· U.E. TS-1 Switch TS-2 .....	P764321-972
· U.E. Chart Drive (60 Cycle) .....	P764321-967
· U.E. Chart Drive (50 Cycle) .....	P764322-096
· U.E. Recorder/Controller .....	P760265-566

## Eagle Small & Medium EO (Gas) Sterilizers:

Existing Bristol Recorder/Controller .....	P141198-272
U.E. Recorder Replacement Kit .....	P764321-966
U.E. Recorder - Available Spare Parts:	
· U.E. Charts (Box of 100) °F .....	P755715-857
· U.E. Charts (Box of 100) °C .....	N/A
· U.E. Pen Arm Asm .....	N/A
· U.E. Pens (Red) .....	P764321-958
· U.E. Pens (Blue) .....	P764321-742
· U.E. Thermal Asm .....	P764321-971
· U.E. TS-1 Switch TS-2 .....	N/A
· U.E. Chart Drive (60 Cycle) .....	P764321-967
· U.E. Chart Drive (50 Cycle) .....	P764322-096
· U.E. Recorder/Controller .....	P141198-272

## Medallion Amscomatic Washer/Sterilizers and Sterilizers:

Existing Bristol Recorder/Controller .....	P760265-564
U.E. Recorder - Available Spare Parts:	
· U.E. Charts (Box of 100) °F .....	P056396-763
· U.E. Charts (Box of 100) °C .....	P056396-764
· U.E. Pen Arm Asm .....	NLA
· U.E. Pens .....	P764321-741
· U.E. Thermal Asm .....	P764321-970
· U.E. TS-1 Switch TS-2 .....	N/A
· U.E. Chart Drive (60 Cycle) .....	P764321-967
· U.E. Chart Drive (50 Cycle) .....	P764322-096
· U.E. Recorder/Controller .....	P760265-564

## REMANUFACTURED STERILIZER DOORS

STERILIZER DOORS		
<u>MODEL/SIZE</u>	<u>REBUILD</u>	<u>NEW</u>
16x16 Medallion Eagle	P764328-576	P760265-258
16x16 Medallion Washer/Sterilizer	NLA	NLA
16x16 Eagle W/S	P764328-578	P760265-268
20x20 Medallion Eagle	P764328-579	P760265-260
24x36 Medallion Power Door	NLA	----
24x36 Eagle Power Door	NLA	----
24x36 Medallion & Eagle (1/2" Thk. Door Plate)	P764324-182	----
24x36 Eagle (5/8" Thk. Door Plate)	P764330-192	P762650-682

### IMPORTANT NOTE:

Door handles and covers are not included in Rebuilding, unless otherwise specified.

**\*\* THIS NOTIFICATION SUPERSEDES ALL PREVIOUS NOTICES. \*\***

### DOOR EXCHANGE PROGRAM

- Enter an order specifying proper part number for rebuilt door.
- Under shipping instructions, enter a note stating "Exchange Program", (if it is an exchange).
- The rebuilt door will be forwarded to the customer in a shipping container designed for return of the existing door.
- If an exchange door is not available, customer's door will have to be rebuilt.

### DOOR REBUILDING PROGRAM

- A Code "D" order and customer purchase order is required to rebuild a customer door.
- Serial number of the sterilizer should be noted on the Code "D" order for proper identification.
- Socket plates, arms and hinges are not included when rebuilding customer door (except at additional cost).
- Door plates are not replaced if they are defective.

### GENERAL NOTES

- B.E.R. (Black End Ring) and AMSCO Cylindrical doors are not included in this program. Small rectangular "L" and "K" plus "MB" doors are no longer available.
- Shipping charges are not included in the pricing.
- Code "C" orders are unacceptable.
- Used sterilizer doors are needed to supplement our rebuilding program. All sizes will be considered and their condition must be acceptable.

## DOOR GASKETS

---

<u>Part Number</u>	<u>Description</u>
P028360-091 .....	24x24" SQUARE
P038965-091 .....	PORTACLAVE
P055676-091 .....	36x42" RECT.
P055677-091 .....	42x48" RECT.
P074366-091 .....	12x16" RECT.
P074367-091 .....	16x16" SQUARE
P074368-091 .....	12" CYLINDRICAL
P074369-091 .....	16" CYLINDRICAL
P074370-091 .....	20" CYLINDRICAL
P074374-091 .....	8x8" SQUARE
P074375-091 .....	11x11" SQUARE
P074376-091 .....	24x24" SQUARE
P074377-091 .....	24x36" RECT.
P074380-091 .....	36x42" RECT.
P074381-091 .....	42x48" RECT.
P074382-091 .....	48x54" RECT.
P074387-091 .....	16" CYLINDRICAL
P074388-091 .....	20" CYLINDRICAL
P074395-091 .....	12" CYLINDRICAL
P074397-091 .....	16" CYLINDRICAL
P074398-091 .....	20" CYLINDRICAL
P078767-091 .....	20x20" SQUARE
P092229-001 .....	28x67" WEDGE
P097487-091 .....	24x36" RECT.
P098051-091 .....	AMSCOMATIC
NLA .....	VERN 20X20" DOOR GASKET
P101007-553 .....	VERN 16X16" CYLINDRICAL DOOR GASKET
P754872-091 .....	20" CYLINDRICAL
P754874-091 .....	601R DYNA
P754875-091 .....	8x8" SQUARE (8816)
P754877-091 .....	24x36" RECT. ME-UME
P754878-091 .....	24x24" SQUARE
P754879-091 .....	24x36" SQUARE
NLA .....	10x16" PORT. ETO
NLA .....	VERN 16" CYLINDRICAL
P115600-001 .....	MEDALIST
P129373-376 .....	CENTURY MEDIUM 26x37.5
P134473-135 .....	STAFF DOOR GASKET
P136800-905 .....	180 STERILIZER
P146653-007 .....	DOOR SEAL (900/901 STERILIZER)
P146660-004 .....	DOOR GASKET (3017)

## DOOR GASKETS

---

<u>Part Number</u>	<u>Description</u>
<b>NLA</b> .....	VERN 24x36
<b>P331008-091</b> .....	62x84" RECT.
<b>P331996-091</b> .....	42x84" RECT.
<b>P426637-261</b> .....	EAGLE 10 STERILIZER
<b>P615851-001</b> .....	30x67" WEDGE
<b>P754871-091</b> .....	16" CYLINDRICAL

*Prices subject to change without notice.*

# VALVE CHART FOR EAGLE 2000 STERILIZERS

	Small Gravity								
	Medium Gravity								
	Small Grav/Iso								
	Med Grav/Iso								
	Small Vac								
	Med Vac								
	Small IWS								
	Small Gas								
	Med Gas								
	Med Vac/Gas								
S-1 Filtered Air					A	A		A	A
S-2 Steam To Ch. (Diaph)	B	C	B	C	B	D	B		D
S-2 Steam To Ch. (Piston)	E	F	E	F	E	G	E		G
S-3 Fast Exhaust	C	C	C	C	C	D	H	C	D
S-4 Exhaust Cooling	I	J	I	J	J	J	I		J
S-5 Steam To Ejector	C	C	C	C	C	C			
S-7 Vacuum Water (Pump)					J	J		I	I
S-8 Vacuum Water (Iso)			K	L					
S-9 Steam To Jacket			E	C					C
S-10 Steam To Ejec Bypass			M	K					
S-11, S-20 Gas Charge								N	N
S-12 Condition. Steam								O	M
S-14 Gas Heater Steam								M	M
S-15 Jacket Drain									C
S-16 Jacket Trap Isolate									B
S-18 Jacket Heat Gas								K	K
S-19 Jacket Cooling									I
S-21 Main Water Supply							P		
IJ-21 Detergent Injector							Q		
S-23 Chamber Drain							R		
S-25 Chamber Water (Fill)							P		
S-26 Steam To Jet Comp							S		
+S-30 Redundant Gas Charge								N	N
+S-31 Gas Slow Exst.								C	C
*Air Vent								T	T
Drain							U		
*S-1 Filtered Air								V	V

+On units shipped after 1/1/87.

\*\*For units with water actuated air vent valves.

## VALVE CHART FOR EAGLE 2000 STERILIZERS (Cont'd.)

	VALVE PART NO.	KIT PART NO.	SIZE NPT
A	P084444-002	P764078-001	3/4
B	P083264-002	P764070-001	3/8
C	P083228-002	P764070-001	3/4
D	P083229-002	P764071-001	1
E	P150822-309	P764317-687	3/8
F	P150822-310	P764317-687	3/4
G	P150822-311	P764317-688	1
H	P083264-003	P764070-001	3/8
I	P083263-002	P764073-001	3/8
J	P083261-002	P764072-001	3/4
K	P084430-002	P764077-001	3/8
L	P093650-002	P764137-001	3/4
M	P083540-002	P764076-001	1/4
N	P150822-934	P764322-880	1/4
O	P129220-001	P764079-001	1/4
P	P083875-002	P764188-001	3/4
Q	P084331-002	N/A	1/2
R	P083475-002	P764075-001	1/4
S	P083874-002	P764187-001	3/8
T*	P093908-334(old)+ P093909-853(new)++	P764320-609 P764323-252	
U	P083885-001	P764189-001	
V**	P056396-480	P764323-939	

*\*Refer to Service Bulletin 89051-I.*

*\*\* If unit has P093900-648, Catalog # J8320B83MO, replace with new valve P056396-480 with stainless-steel seat.*

*+Old style valve, Catalog HV226-187-2*

*++New style valve, Catalog HV238-174-1*

**NOTE: SERVICE BULLETIN #87036-I  
VARISTOR P150778-001**

## VALVE CHART FOR EAGLE 3000 STAGE 2 STERILIZERS

	Small Gravity				
	Medium Gravity				
	Small Vac				
	Medium Vac				
	Small IWS				
S-1 Filtered Air			A	A	
S-2 Steam To Ch. (Piston)	E	G	E	G	E
S-3 Fast Exhaust	C	C	C	D	B
S-4 Exhaust Cooling	I	J	J	J	I
S-5 Steam To Ejector	C	C	C	C	
S-7 Vacuum Water (Pump)			J	J	
S-9 Steam To Jacket	E	C	B	C	
S-11, S-20 Gas Charge					
S-12 Condition. Steam					
S-14 Gas Heater Steam					
S-18 Jacket Heat Gas					
S-21 Main Water Supply					I
IJ-21 Detergent Injector					Q
S-23 Chamber Drain					I
S-25 Chamber Water (Fill)					P
S-26 Steam To Jet Comp					K
S-30 Redundant Gas Charge					
S-31 Gas Slow Exst.					
S-35 Purge Valve					
S-41 Exhaust Valve (Disposer Only)					
S-42 Fast Exhaust (Disposer Only)					
S-50 Jacket Drain					
CV-21 Chamber Drain					S
Air Vent Valve					
S9W (Auto Utility Option)	J	J	J	J	



**VALVE CHART FOR EAGLE 3000 STAGE 2 STERILIZERS (Cont'd.)**

	VALVE PART NO.	P-KIT PART NO.	SIZE NPT
A	P084444-002	P764078-001	3/4
B	P083264-002	P764070-001	3/8
C	P083228-002	P764070-001	3/4
D	P083229-002	P764071-001	1
E	P150822-309	P764317-687	3/8
G	P150822-311	P764317-688	1
I	P083263-002	P764073-001	3/8
J	P083261-002	P764072-001	3/4
K	P084430-002	P764077-001	3/8
P	P083875-002	P764188-001	3/4
Q	P084331-002	N/A	1/2
S	P083885-001	P764189-001	1-1/4

## VALVE CHART FOR RENAISSANCE STERILIZERS

	Small Gravity			
	Medium Gravity			
	Small Vac			
	Medium Vac			
S-1 Filtered Air			A	A
S-2 Steam To Ch. (Piston)	E	F	E	F
S-3 Fast Exhaust	C	C	C	D
S-4 Exhaust Cooling		B	B	
S-5 Steam To Ejector		C		C
S-7 Vacuum Water	B		B	B
S-9 Steam To Jacket	E	F	E	F

	VALVE PART NO.	P-KIT PART NO.	SIZE NPT
A	P093911-329	P764324-895	3/4
B	P093911-328	P764072-001	3/4
C	P093910-477	P764070-001	3/4
D	P093911-327	P764071-001	1
E	P093910-650	P764317-687	3/8
F	P093911-331	P764317-688	1

## VALVE CHART FOR EAGLE 3000 STAGE 3 STERILIZERS

	Small Gravity									
	Medium Gravity									
	Small Vac									
	Medium Vac									
	Small IWS									
	Small Gas									
	Medium Gas									
	Small Isothermal									
	Med. Isothermal									
S-1 Filtered Air			A	A		F	F			
S-2 Steam To Ch. (Piston)	E	G	E	G	E			E	G	
S-3 Fast Exhaust	C	C	C	D	H	C	C	C	C	
S-4 Exhaust Cooling		J		J	I				J	
S-5 Steam To Ejector	C	C	C	C				C	C	
S-7 Vacuum Water	J		J	J		I	I	J		
S-8 Steam Ejector								E	S	
S-9 Steam To Jacket	E	G	E	G				E	G	
S-10 Ejector Bypass								R	E	
S-11, S-20 Gas Charge						N	N			
S-12 Condition. Steam						O	R			
S-14 Gas Heater Steam						R	R			
S-16 Main - Water	D	D	D	D				D	D	
S-18 Jacket Heat Gas						E	E			
S-21 Main Water Supply					I					
IJ-21 Detergent Injector					Q					
IJ-22 Injector					K					
S-23 Chamber Drain					I					
S-25 Chamber Water (Fill)					P					
S-26 Steam To Jet Comp					E					
S-27 Multiport Shutoff	B	B	B	B				B	B	
S-30 Redundant Gas Charge						N	N			
S-31 Gas Slow Exst.							C			
S-35 Purge Valve						N	N			
S-41 Exhaust Valve (Disposer Only)						C	C			
S-42 Fast Exhaust (Disposer Only)							C			
S-50 Jacket Drain						C	C			
V-23 Chamber Drain					M					
Air Vent Valve						L	L			

*NOTE: S13 Valve used for Optimal Solution Cooling, refer to next page.*

## VALVE CHART FOR EAGLE 3000 STAGE 3 STERILIZERS (Cont'd.)

	VALVE PART NO.	P-KIT PART NO.	SIZE NPT
A	P093911-329	P764324-895	3/4
B	P093911-364	P764324-090	1/8
C	P093910-477	P764070-001	3/4
D	P093911-327	P764071-001	1
E**	P093910-650	P764317-687	3/8
F*	P093910-476	P764323-939	1/8
G	P093911-331	P764317-688	1
H	P093911-332	P764070-001	3/8
I	P093910-550	P764073-001	3/8
J	P093911-328	P764072-001	3/4
K	P084318-001	N/A	1/8
L	P093909-853	P764323-252	3/4
M	P083885-001	P764189-001	1-1/4
N	P093910-481	P764322-880	1/4
O	P093910-478	P764079-001	1/4
P	P093911-333	P764188-001	3/4
Q	P084331-002	N/A	1/2
R	P093910-479	P764076-001	1/4
S	P093911-551	P764317-687	3/4

\* Units shipped before 6/90 used a brass type valve, catalog #SC8320G83. Replace with new valve P093910-476 (stainless-steel type).

\*\* Units shipped before 8/90 used on P093910-798, catalog HV258-902-1. Replace with new valve P093910-650.

S13 used on Optimal Cooling units:  
P093910-478      1/4" Valve  
P764079-001      Valve Kit  
P764323-942      Coil

**VALVE CHART FOR EAGLE 2400 AND 3400  
AMSCOMATIC WASHER/STERILIZERS**

		<b>2400</b>	<b>3400</b>
V-1	MAIN WATER FILL	C	I
V-2	TAP WATER FILL	C	K
V-3	WASH PUMP WATER	A	G
V-4	DETERGENT	B	H
V-5	PURE WATER FILL	C	K
V-6	STEAM TO SUMP	D	J
V-7	DRAIN	F	L
V-8	STEAM TO CHAMBER	E	K
V-9	WATER TO EJECTOR	D	J
V-10	VACUUM EJECTOR	D	J
V-11	VACUUM BREAK	D	J
V-12	STEAM TRAP		J
V-13	WATER TO CONDENSER	A	
	PNEUMATIC VALVE	M	
	PNEUMATIC VALVE		N

*\* NOTE: This repair kit is for use on actuator cat. no. A412 or 412 SRVO only and is NLA. It will be necessary to replace the entire actuator with the newer style B412 series actuator (P764326-717). A actuator mounting kit will also be required. To determine the correct mounting kit to be used, both the pipe thread size (NPT) and the model/series of the valve the actuator is to be used on is required. For 1/2 NPT Series 76 or 1/4 NPT and 3/8 NPT Series 86, use mounting kit P764326-714. For 3/4 NPT Series 76 or 1/2 NPT Series 86, use mounting kit P764326-715. For 1 NPT Series 76 or 3/4 NPT Series 86 or 1/2 NPT Series 96, use mounting kit P764326-716. Repair kit for actuator B412 (P764326-717), use P764329-295.*

# VALVE CHART FOR EAGLE 2400 AND 3400 AMSCOMATIC WASHERS/STERILIZERS (Continued)

	SIZE NPT	VALVE	VALVE REPAIR KIT	VALVE BODY COMPLETE	ACTUATOR	ACTUATOR REPAIR KIT	ACTUATOR GASKET KIT
A	1/4	P041626-091 (TYPE J-218) (TYPE GP-218)	NLA P764319-924				
B	3/8	P056396-203	P764319-127				
C	1	P080028-091 (TYPE 810P) (TYPE RB41P8)	P759494-001 P764319-946				
D	1/2	P150822-928	P764317-587	P764317-863	P764326-717	See NOTE*	NLA
E	3/4	P150822-929	P764317-588	P764316-938	P764326-717	See NOTE*	NLA
F	1-1/4	P056399-117	P764317-590	P764317-043	NLA	P764317-579	P764317-579
G	1/4	P093909-883	P764323-833				
H	3/8	P093910-599	P764323-834				
I	1	P093910-600	P764323-835				
J	1/2	P150822-928	P764317-587	P764317-863	P764326-717	See NOTE*	NLA
K	3/4	P150822-929	P764317-588	P764316-938	P764326-717	See NOTE*	NLA
L	1-1/4	P056399-117	P764317-590	P764317-043	P764326-854	P764317-579	P764317-579
M	1/8	P764317-658	P764324-745				
N	1/8	P093910-591					

NOTE: See previous page for note on repair kit.

## SYSTEM 2S/LAB 90 SOLENOID VALVES

<b><u>Piping</u></b>	<b><u>Part Number</u></b>	<b><u>Description</u></b>
S1	P093922-054	Valve, Sol. 1/4" NPT (Air Break)
	P764328-740	Kit, Repair
	P764328-737	Coil
S2	P093922-033	Valve, Sol. 1/4" NPT (Chamber Strm.)
	P764328-732	Kit, Repair
	P764328-733	Coil
S3	P093922-032	Valve, Sol. 1/2" NPT (Chamber Exh.)
	P764328-738	Kit, Repair
	P764328-739	Coil
S7	P093922-031	Valve, Sol. 1/2" NPT (Cooling Wtr./Ejector)
	P764328-736	Kit, Repair
	P764328-737	Coil
S35	P093922-033	Valve, Sol. 1/4" NPT (Seal Strm.)
	P764328-732	Kit, Repair
	P764328-733	Coil
S37	P093922-033	Valve, Sol. 1/4" NPT (Seal Exhaust)
	P764328-732	Kit, Repair
	P764328-733	Coil
S50	P093922-031	Valve, Sol. 1/2" NPT (Filtered Water)
	P764328-736	Kit, Repair
	P764328-737	Coil

## CENTURY/MILLENNIUM STERILIZERS VALVES

<u>Piping</u>	<u>Part Number</u>	<u>Description</u>
<b>Steam Supply Manifold (P146657-787)</b>		
	P764326-475	Manifold Repair Kit (less seats)
	P764326-476	Seat Package (complete)
S1	N/A	Valve, Solenoid (Air Break)
	P764326-477	Kit, Repair
	P764326-632	Coil
	NLA	Seat
S2, S9	N/A	Valve, Solenoid (Steam to Chamber/ Steam to Jacket)
	P764326-478	Kit, Repair
	P764323-941	Coil
	NLA	Seat
S35	N/A	Valve, Solenoid (O.E. Door Seal)
	P764326-479	Kit, Repair
	P764323-941	Coil
	NLA	Seat
<b>Exhaust Manifold (P146657-788)</b>		
	P764326-480	Manifold Repair Kit (less seats)
	P764326-481	Seat Package (complete)
S3	N/A	Valve, Solenoid (Fast Exhaust)
	P764326-482	Kit, Repair
	P764323-941	Coil
	NLA	Seat
S37	N/A	Valve, Solenoid (O.E. Door Exhaust)
	P764326-479	Kit, Repair
	P764323-941	Coil
	NLA	Seat
S40	N/A	Valve, Solenoid
	P764326-483	Kit, Repair
	P764323-941	Coil
	NLA	Seat
TR1	P764315-209	Steam Trap, Chamber (cap only)
	P764080-001	Kit, Repair
CK4,6,FC2		Check Valve/Flow Control Valve
	P764326-485	Kit, Repair
	NLA	Seat (CK4, CK6)
CK8		Check Valve
	P764326-484	Kit, Repair
	NLA	Seat



## CENTURY/MILLENNIUM STERILIZERS VALVES (continued)

<u>Piping</u>	<u>Part Number</u>	<u>Description</u>
<b>Water Supply Manifold (P146657-789)</b>		
	P764326-619	Manifold Repair Kit (less seats)
	P764326-620	Seat Package (complete)
S4, S7	N/A	Valve, Solenoid (Cooling Water Vacuum Water)
	P764326-621	Kit, Repair
	P764323-940	Coil
	NLA	Seat
CK5	N/A	Check Valve
	P764326-622	Kit, Repair
	NLA	Seat
<b>Double Door Option Manifold (P136809-791)</b>		
	P764326-486	Manifold Repair Kit (less seats)
	P764326-487	Seat Package (complete)
S38	N/A	Valve, Solenoid (N.O.E. Door Exhaust)
	P764326-479	Kit, Repair
	P764323-941	Coil
	NLA	Seat
FC3, CK7	N/A	Flow Control Valve/Check Valve
	P764326-485	Kit, Repair
	NLA	Seat (CK7)

**AMSCO CENTURY MEDIUM STERILIZER SOLENOID VALVES  
(26x37.5 CHAMBER, PREVACUUM & SFPP MODELS)**

<b><u>Piping</u></b>	<b><u>Part Number</u></b>	<b><u>Description</u></b>
S1 (NO)	P093911-329 P764324-895 P764323-941	Sol. Valve, Filtered Air, 3/4"NPT Kit, Repair Coil
S2 (NC)	P093911-331 P764317-688 P764323-941	Sol. Valve, Steam to Chamber, 1"NPT Kit, Repair Coil
S3 (NC)	P093911-327 P764071-001 P764324-600	Sol. Valve, Fast Exhaust, 1"NPT Kit, Repair Coil
S4 (NC)	P093911-328 P764072-001 P764323-940	Sol. Valve, Exhaust Cooling, 3/4"NPT Kit, Repair Coil
S7 (NC)	P093910-479 P764076-001 P764323-942	Sol. Valve, Vacuum Water, 1/4"NPT Kit, Repair Coil
S9 (NC)	P093911-331 P764317-688 P764323-941	Sol. Valve, Steam to Jacket, 1"NPT Kit, Repair Coil
S35 (NC)	P093910-479 P764076-001 P764323-942	Sol. Valve, Steam to Seal (OE), 1/4"NPT Kit, Repair Coil
S36 (NC)	P093910-479 P764076-001 P764323-942	Sol. Valve, Steam to Seal (NOE), 1/4"NPT Kit, Repair Coil
S37 (NC)	P093910-479 P764076-001 P764323-942	Sol. Valve, Seal Retraction (OE), 1/4"NPT Kit, Repair Coil
S38 (NC)	P093910-479  P764076-001 P764323-942	Sol. Valve, Seal Retraction (NOE), 1/4"NPT Kit, Repair Coil
S40 (NO)	P093910-478  P764079-001 P764323-942	Sol. Valve, Chamber Constant Bleed, 1/4"NPT Kit, Repair Coil
S43 (NC)	P093910-479 P764076-001 P764323-942	Sol. Valve, Vacuum Pump Drain, 1/4"NPT Kit, Repair Coil
S46 (NC)	P093910-479 P764076-001 P764323-942	Sol. Valve, Waste Cooling Water, 1/4"NPT Kit, Repair Coil

AMSCO SMALL AND MEDIUM STERILIZER  
**SAFETY VALVES**

TYPE OF UNIT	MEDALLION	EAGLE 2000	EAGLE 3000 STAGE 2	EAGLE 3000 STAGE 3	RENAISSANCE	SMALL CENTURY 116/120	MEDIUM CENTURY 26x37.5	SYSTEM 2S/ LAB 90
SMALL GENERAL PURPOSE	STEAM HEAT P764316-929 ELECTRIC HEAT P150526-001	P150828-473	P150828-473	P150828-473	P150828-473	P093918-013		
SMALL ISOTHERMAL	P764316-929	P150828-473		P150828-473		P093918-013		
SMALL VACAMATIC	P150828-476	P150828-473	P150828-473	P150828-473	P150828-473	P093918-013		P093918-013
SMALL WASHER	P764316-929	P084342-001	P084342-001	P084342-001				
SMALL ETO (CHAMBER)	P150828-475	P150828-475		P150828-475				
SMALL ETO (JACKET)				P150828-474				
MEDIUM GENERAL PURPOSE	P150828-476	P150828-476	P150828-476	P150828-478	P150828-478			
MEDIUM ISOTHERMAL	P150828-476	P150828-476		P150828-478				
MEDIUM VACAMATIC	P150828-476	P150828-476	P150828-476	P150828-478	P150828-478		P093921-266	
MEDIUM ETO (CHAMBER)	P764324-576	P764324-576		P150828-475				
MEDIUM ETO (JACKET)	P150828-476	P150828-473		P150828-474				

FLS safety valve P150828-476

FLJ safety valve P150828-478

## **SAFETY VALVE SPECIFICATIONS**

<b>PART NUMBER</b>	<b>DESCRIPTION</b>
P084342-001	36 PSIG SET PRESSURE, 3/4" NPT, 1" NPT OUTLET
P150526-001	33 PSIG SET PRESSURE, 1/2" NPT, 3/4" NPT OUTLET
P150828-473	40 PSIG SET PRESSURE, 1/2" NPT, 3/4" NPT OUTLET
P150828-474	15 PSIG SET PRESSURE, 1/2" NPT, 3/4" NPT OUTLET
P150828-475	15 PSIG SET PRESSURE, 1/2" NPT, 3/4" NPT OUTLET
P150828-476	40 PSIG SET PRESSURE, 3/4" NPT, 1" NPT OUTLET
P150828-478	40 PSIG SET PRESSURE, 1" NPT, 1-1/4" NPT OUTLET
P764316-929	36 PSIG SET PRESSURE, 1/2" NPT, 3/4" NPT OUTLET
P764324-576	15 PSIG SET PRESSURE, 1/2" NPT, 3/4" NPT OUTLET
P093918-013	50 PSIG SET PRESSURE, 1/2" NPT, 3/4" NPT OUTLET
P150828-476	40 PSIG SET PRESSURE, 3/4" NPT, 1" NPT OUTLET
P150828-478	40 PSIG SET PRESSURE, 1" NPT, 1-1/4" NPT OUTLET
P093921-266	45 PSIG SET PRESSURE, 1" NPT, 1-1/4" NPT OUTLET

**EO EAGLE MONITOR STAINLESS STEEL FITTING  
& COMPONENTS**

FITTING & SIZE	FRONT FERRULE	BACK FERRULE	NUT
P755716-835 Bulkhead Union, 1/8"	P755716-844 Ferrule Front, 1/8"	P755716-843 Ferrule Back, 1/8"	P755716-845 Nut, 1/8"
P755716-836 Bulkhead Union, 1/4"	P200000-208 Ferrule Front, 1/4"	P200000-208 Ferrule Back, 1/4"	P764321-118 Nut, 1/4"
P755716-837 Tee Union, 1/8"	P755716-844 Ferrule Front, 1/8"	P755716-843 Ferrule Back, 1/8"	P755716-845 Nut, 1/8"
P755716-838 Union Elbow, 1/8"	P755716-844 Ferrule Front, 1/8"	P755716-843 Ferrule Back, 1/8"	P755716-845 Nut, 1/8"
P755716-839 Union, 1/8"	P755716-844 Ferrule Front, 1/8"	P755716-843 Ferrule Back, 1/8"	P755716-845 Nut, 1/8"
P755716-840 Reducer, Union 1/8" x 1/4"	P755716-844 Ferrule Front, 1/8"	P200000-208 Ferrule Back, 1/4"	P755716-845 Nut, 1/8" P764321-118 Nut, 1/4"
P755716-841 Male Elbow, 1/8" x 1/4"	P755716-844 Ferrule Front, 1/8"	P200000-208 Ferrule Back, 1/4"	P755716-845 Nut, 1/8" P764321-118 Nut, 1/4"
P755716-842 Male Elbow, 1/4" x 1/4"	P200000-208 Ferrule Front, 1/4"	P200000-208 Ferrule Back, 1/4"	P764321-118 Nut, 1/4"
P755716-847 1/8 x 1/4" Male Connector	P200000-208 Ferrule Front, 1/4"	P200000-208 Ferrule Back, 1/4"	P764321-118 Nut, 1/4"

## SURGICAL PRODUCTS CONSUMABLES

### LIGHTS (ROTOFLEX/CHALLENGE 22/POLARIS)

TRACKS:	P751446-091	GREASE - DOW CORNING
CARRIAGES:	R006400-703	ALVANIA EP-2 GREASE - 16 OZ TUBE
ARM & LIGHTHEADS:	R006400-703	ALVANIA EP-2 GREASE - 16 OZ TUBE
	R005300-542	LOCTITE 609-31 BEARING RETAIN - 50CC BOTTLE
	R005300-548	LOCTITE 271-21 ADHESIVE - 10CC BOTTLE
	R006200-402	LINT FREE WIPES - BOX OF 100
CENTRA 360:	R005300-542	LOCTITE 609-31 BEARING RETAIN - 50CC BOTTLE
	R006400-180	SNOW WHITE PETROLEUM JELLY - 20 OZ CONTAINER

### POLARIS SL LIGHT

ARM & LIGHTHEADS:	R005300-162	RTV 106 ADHESIVE - 3 OZ TUBE
	R005300-542	LOCTITE 609-31 BEARING RETAIN - 50CC BOTTLE
	R005300-545	LOCTITE SEALANT # 222-31 - 50CC BOTTLE
	R005300-548	LOCTITE 271-21 ADHESIVE - 10CC BOTTLE
	R005300-554	LOCTITE 242 - 50ML BOTTLE
	R006200-402	LINT FREE WIPES - BOX OF 100
	R006400-703	ALVANIA EP-2 GREASE
	R006400-826	MO-LITH # 2 LUBRIPLATE - 14-1/2 OZ TUBE
	R006400-880	MOLYKOTE # 321-R - 13 OZ SPRAY CAN
CENTRA 360:	P0431215-091	PLIOBOND
	R005300-542	LOCTITE 609-31 BEARING RETAIN - 50CC BOTTLE
	R005300-545	LOCTITE SEALANT # 222-31 - 50CC BOTTLE

## **SURGICAL PRODUCTS CONSUMABLES (Cont'd.)**

### POLARIS SL LIGHT

CENTRA 360: (Cont'd.)	R006400-180	SNOW WHITE PETROLEUM JELLY - 20 OZ CONTAINER
--------------------------	-------------	---

### EXAMINER 10 LIGHT

BRACKET ASSEMBLY:	R006400-703	ALVANIA EP-2 GREASE - 16 OZ TUBE
----------------------	-------------	-------------------------------------

SUSPENSION ARM:	R006400-310	ORE-LUBE NEPTUNE 7 GREASE - 14-1/2 OZ TUBE
--------------------	-------------	---

### 900 LIGHT

	R006400-310	ORE-LUBE NEPTUNE 7 GREASE - 14-1/2 OZ TUBE
--	-------------	---

### 180 LIGHT

LIGHTHEAD ASSEMBLY	R005300-548	LOCTITE 271-21 ADHESIVE - 10CC BOTTLE
-----------------------	-------------	--

	R005300-554	LOCTITE 242 - 50ML BOTTLE
--	-------------	---------------------------

ARM AND LIGHTHEAD	R005300-563	RTV 108 GE ADHESIVE - 3 OZ TUBE
----------------------	-------------	------------------------------------

	R006400-703	ALVANIA EP-2 GREASE - 16 OZ TUBE
--	-------------	-------------------------------------

	R006400-853	JET-LUBE AP-1W WHITE GREASE - 1 PINT CAN
--	-------------	---

### GEMINI LIGHT

	R005300-545	LOCTITE SEALANT # 222-31 - 50CC BOTTLE
--	-------------	---

	R005300-563	RTV-108 GE ADHESIVE - 3 OZ TUBE
--	-------------	------------------------------------

	R006400-703	ALVANIA EP-2 GREASE - 16 OZ TUBE
--	-------------	-------------------------------------

	R006400-853	JET-LUBE AP-1W WHITE GREASE - 1 PINT CAN
--	-------------	---

	R005300-554	LOCTITE 242 - 50ML BOTTLE
--	-------------	---------------------------

	R006400-880	MOLYKOTE # 321-R - 13 OZ SPRAY CAN
--	-------------	---------------------------------------

## **SURGICAL PRODUCTS CONSUMABLES (Cont'd.)**

### **QUANTUM LIGHT**

ARM AND LIGHTHEAD:	R005300-286	ORE-LUBE ANTI-SEIZE COM- POUND - 14-1/2 OZ CARTRIDGE
	R005300-545	LOCTITE SEALANT # 222-31 - 50CC BOTTLE
	R005300-557	LOCTITE # 495 ADHESIVE - 1 OZ CONTAINER
	R005300-890	LOCTITE # 262 THREAD LOCKER - 50ML BOTTLE
	R006400-180	SNOW WHITE PETROLEUM JELLY - 20 OZ CONTAINER
CENTRA 360:	R006400-703	ALVANIA EP-2 GREASE - 16 OZ TUBE
	R005300-542	LOCTITE 609-31 BEARING RETAIN - 50CC BOTTLE
	R005300-545	LOCTITE SEALANT # 222-31 - 50CC BOTTLE
VARIAC CONTROL:	R006400-180	SNOW WHITE PETROLEUM JELLY - 20 OZ CONTAINER
	R005300-545	LOCTITE SEALANT # 222-31 - 50CC BOTTLE

### **TABLES**

2080 SERIES:	R006400-826	MOLY-LUBRIPLATE TYPE MSHD #2 - 5 LB CAN
	P754811-091	LUBRIPLATE TYPE 630AA - 5 LB CAN
3080 SERIES:	R005300-286	ORE-LUBE ANTI-SEIZE COM- POUND - 14-1/2 OZ CARTRIDGE
	R005300-542	LOCTITE 609-31 BEARING RETAIN - 50CC BOTTLE
	R005300-545	LOCTITE SEALANT # 222-31 - 50CC BOTTLE
	R005300-554	LOCTITE 242 - 50ML BOTTLE
	R005300-563	RTV 108 GE ADHESIVE - 3 OZ TUBE



## **SURGICAL PRODUCTS CONSUMABLES (Cont'd.)**

### TABLES (Cont'd.)

3080 SERIES: (Cont'd.)	P764326-267	680 CHEVRON GEAR LUBE (4 OZ)
	R005300-587	ADHESIVE ACCELERATOR VERSILOK - 1 GALLON
	R006200-402	LINT FREE WIPES - BOX OF 100
	R006400-826	MO-LITH # 2 LUBRIPLATE - 14-1/2 OZ TUBE

### HYDRAULIC OIL FOR TABLES - TWO 2 QT. CONTAINERS

#### 3080RC & 3080RL:

P764322-636	REFINED BY MONTGOMERY
-------------	-----------------------

#### 2080M, 2080L, 2080RC:

P764322-636	REFINED BY MONTGOMERY
-------------	-----------------------

1080M:	P764257-002	SOLNUS # 750
--------	-------------	--------------

1080L:	NLA	SUNOCO SOLNUS # 300
--------	-----	---------------------

#### ORTHOGRAPHIC:

P082182-001	SHELL #X100SAE 40
-------------	-------------------

#### ORTHOGRAPHIC 2 & ORTHOVISION:

P764322-636	REFINED BY MONTGOMERY
-------------	-----------------------

#### DELIVERY 4:

P764322-636	REFINED BY MONTGOMERY
-------------	-----------------------

#### ALBEE FRACTURE TABLE:

P764257-002	SOLNUS # 750
-------------	--------------

GU110:	P764322-636	REFINED BY MONTGOMERY
--------	-------------	-----------------------

#### 500 O.B. TABLE:

P764257-002	SOLNUS # 750
-------------	--------------

#### 800 & 800W (MANUAL):

P764257-002	SOLNUS # 750
-------------	--------------

#### 800 & 800W (ELECTRIC):

NLA	SUNOCO SOLNUS # 300
-----	---------------------

## TOOL LIST

<u>PART NUMBER</u>	<u>DESCRIPTION</u>
P054145-045	Wrench, Cylinder Adjusting
P752608-045	Wrench, For Locknut (Sterilizer Doors)
P752919-091	Wrench, Monkey - General Duty 11" Long, Open 2-3/4"
P752967-091	Wrench, Pipe - Heavy Duty 6"
P752968-091	Ratchet Kit, Midget
P752969-091	Hooks, Flexible Packing - 3 x 5" Long
P752970-091	Kit, Spline Key
P752973-091	Scraper, Bearing 3-Cornered
P752974-091	Fuse, Puller
P752975-091	Meter, Volt/Amp. - Amprobe RS-300 Carrying Case, Instruction Manual & Test Leads Incl.
P752980-091	Wrench, Thin Pattern - 2" Jaw
P752987-091	Adjusting Tool - Cat. No. 146 (SAGA)
P754219-091	Tag, Safety
P754763-091	Wrench, Spanner (1080 Table)
P755523-001	Wrench, Pedestal Adjusting Screw
P755715-123	Kit, Cleaning
P755715-905	*Multimeter (Fluke 87)
P755716-176	Tool Insertion/Extraction
P757264-045	Wrench, Spanner (2080 Table)
P757290-091	Scale, Spring - 0-50 Lbs. Capacity
P758229-001	Tool, Bearing Insertion for 2080 Tables
P759729-001	Wrench, Torque - 3/8" 0-300 In./Lbs.
P759730-001	Wrench, Torque - 1/2" 0-150 Ft./Lbs.
P764330-525	*Kit, Light Metering
P764317-559	Case, Carrying
P764319-361	*Gauge, Pressure W/Case
P764328-379	Tool Box - 18x16x7" (Dk. Grey - New) 2 per Tech
P764328-381	• Handle, Replacement
P764328-382	• Latch, Replacement
P755716-396	• Pallets, Replacement
NLA	Tool Box - 18x14x9-1/2" (Grey/Silver - Old) 2 per Tech
P764328-381	• Handle, Replacement
P764328-383	• Latch, Replacement
P755716-275	• Pallets, Replacement
P764322-660	Tool Kit (3080 Table)
NLA	*Pyrometer w/Case
NLA	*Pyrometer, Hand Held
P764327-796	Cart, Tool Case
P764330-629	• Top Handle
P764330-628	• Wheel, 6"
P764322-952	Wrench, Special Allen
P764326-076	Hand Held Viewer
P764328-215	• Lamp, 2.5V
P764328-216	• Power Cord

**\*NOTE:** These tools need calibration on a six month basis. Advise calibration control number.

## TOOL LIST (continued)

<u>PART NUMBER</u>	<u>DESCRIPTION</u>
P764328-167	Kit, Calibration "Tech Only"
P000939-091	• Red Bsg, 1/4NPT, 1/8NPT (qty 2)
P028899-091	• Nipple, 1/8 x 1", Brass (qty 2)
P129353-531	• Male Connector, 1/4ODT x 1/8NPT (qty 2)
P764317-559	• Carrying Case
P764326-254	• Hose Assembly
P764327-984	• Heise Pte Base Unit
P764327-985	• Pressure Sensor, 100 psi
P764327-986	• Pressure Sensor, 60 psi
P764328-805	• Pressure Sensor Module
P764329-246	• Snubber Fitting
R003500-735	• Nylon Tubing, 1/4OD x 180ID (qty 3)
P764321-342	Plug, Temp. Probe (Printcon Test Plug)

## MISCELLANEOUS TOOLS

P763645-001	Helmet, Safety
NLA	Board, Clip
NLA	Panel, Microfiche

## EO/SAFETY RELATED TOOLS

P764316-329	Headgear
P764316-330	Shield
P764316-331	Gas Mask
	Replacement Face Pieces:
P755716-660	• Small
P755716-661	• Medium
P755716-662	• Large
P764316-332	Canister
P764316-690	Jumpsuit (Medium)
P764316-691	Jumpsuit (Large)
P764316-692	Jumpsuit (Extra Large)
P764325-929	Jumpsuit (Extra-Extra Large)
P764328-797	EO Gas Detector, "Tech Only"
NLA	Lens Kit
NLA	Cleaner-Sanitizer 25 1-oz. Pkgs.
P764323-796	Gloves, Best Butyl Rubber (Pair)
NX-001	EO Badge
P150828-230	Mitt (Oven)
P764324-568	Sleeve, Protective Kit (Kevlar)
NLA	EO Safety Kit
P764327-651	Ear Muffs
P764321-583	Wrist Strap
P764329-984	Knee Pads (Pr.)

## SERVICE LITERATURE

DESCRIPTION	PART NUMBER
Service Kit Listing Book	P764322-055
Section "Q" Valves & Traps	P764322-401
Catalog, Table Accessory Components	P764322-008
Service Bulletin Book	SV-6100
Serviceman's Tool Record	2-40-9
Service Quick Reference Handbook	SV-5042
Microfiche Set, "Tech Only"	NLA

## STERIS TECHNICIAN LABELS, TAGS AND SIGNS

DESCRIPTION	PART NUMBER
Service Vehicle Sign	P755715-579
Preventive Maintenance Label	P764329-231
Calibration Label	NLA
Decontamination Tag	P764327-797
Return for Evaluation Label	P764329-779
Montgomery Core Return Label	P764327-559
"This Equipment is Not Working Now..." Tag	SV-5099
"This Equipment Ready to Work..." Tag	SV-5100
"Danger, Ethylene Oxide..." Sign	P755715-467
"Danger, Potential Cancer..." Sign	P755715-859
"Notice, No Admittance..." Sign	P755715-465
"Warning, Sterilization Maintenance..." Sign	P755715-851
"Caution, Envirogard - Room Exhaust..." Sign	P755715-466
"ETO, Emergency Locker" Sign	NLA

# PRINTER PAPER MATRIX

Product	Type of Printer		Printer Paper		Ribbon Cartridge	Printer
	Thermal	Impact	1-Ply	2-Ply		
Eagle 2000 - Olivetti		X	P129356-116	P129356-117	P129356-118	N/A
Eagle 2000 - Epson		X	P764324-057	P150828-420	P150828-440	P056401-143
Eagle 3000 - Stage 2	X		P129359-008	N/A	N/A	P764324-604
Eagle 3000 - Stage 3	X		P129359-008	N/A	N/A	P755716-005
Eagle 3000 - Stage 3 *		X	P764324-057	P150828-420	P150828-440	P093918-051
Eagle 3000 - Stage 3 **		X	P764324-057	P150828-420	P150828-440	P338517-173
Reliance 3400X		P129359-008	N/A	N/A	P755716-005	
Eagle Ten	X		P129359-008	N/A	N/A	P755716-005
Reliance Washers	X		P129359-008	N/A	N/A	P755716-005
Reliance Washers **		X	P764324-057	P150828-420	P150828-440	P338517-173
Basil Washers X		P129359-008	N/A	N/A	P755716-005	
Millennium/Century Series		X	P129362-819	P150828-420	P150828-440	P093918-051
Eagle Model 3017		X	P129362-819	N/A	P150828-440	P093911-502
Eagle 3000SL	X	P129362-819	N/A	P150828-440	P093911-502	
VHP 1000	X		P129359-008	N/A	N/A	P755716-005
VHP 1000		X	P093914-204	P150828-420	P150828-440	P093914-003
VHP 100		X	P093914-204	P150828-420	P150828-440	P093914-003
DV 1000		X	P093914-204	P150828-420	P150828-440	P093914-003
DV 30		X	P093914-204	P150828-420	P150828-440	P093914-003
System 1, Models 89, 90	X		P129359-008	N/A		
System 1, Models 99		X	P129362-819	N/A	P150828-440	
System 2S/Lab 90		X	P129362-819	N/A	P150828-440	

\* Impact printer without cable, 15 position header (for units built after 1/96) - 8 dip switches on P/C board

\*\* Impact printer with cable attached, 15 pin connector - 4 dip switches on P/C board

**SMALL AND MEDIUM STERILIZERS  
VACUUM BREAKERS/REPAIR KITS**

<b>Product</b>	<b>Medallion Sterilizer</b>	<b>Eagle 2000 Sterilizer</b>	<b>Eagle 3000 Stage 2 Sterilizer</b>	<b>Eagle 3000 Stage 3 Sterilizer</b>	<b>Renaissance Sterilizer</b>
Gravity - Small	1.	1.	4.	4.	4.
Gravity - Medium	1.	4.	4.	4.	4.
Vacumatic - Small	3.	2.	4.	4.	4.
Vacumatic - Medium	1.	2.	4.	4.	4.
Ethylene Oxide - Small	1.	4.		2.	
Ethylene Oxide - Medium	1.	1.		2.	
Washer/Sterilizer	1.	4.	4.	4.	
Lab/Isothermal - Small	1.	1.		4.	
Lab/Isothermal - Medium	1.	4.		4.	

1. P090233-091 3/8" NPT
  - P752477-091
2. P077021-091 1/2" NPT\*
  - P757815-091(V370A) Sloan
  - P764323-984 (Zum) N.L.A.
  - P764326-720 (38-103) Conbraco
  - P764329-949 (38-203) Conbraco
3. P077022-091 3/8" NPT\*
  - P751077-091 (V370A) Sloan
  - P764326-202 (38-102) Conbraco
  - P764329-948 (38-202) Conbraco
4. P077023-001 3/4" NPT\*
  - P752735-091 (370A) Sloan
  - P764323-804 (Zum) N.L.A.
  - P764326-721 (38-104) Conbraco
  - P764329-950 (38-204) Conbraco

\* Determine manufacture model number of vacuum breaker to order proper repair kit.

# ELECTRO-STEAM STEAM GENERATOR HEATING ELEMENTS

Size/Type	Mat'l	Voltage						Indirect (Stm to Stm)	Heater Gasket
		208	240	380	415	480	575/600		
15AS Non-Solid State*	C/S	P752709-091	P077318-091	P756204-091	P759246-001	P754749-091	P764322-469		P764321-472 (Set of 3)
	S/S	P764315-577	P764315-578	P764315-579	P764315-579	P764315-580		P764321-469	
20AS Non-Solid State*	C/S	P078826-091	P075425-091	P764317-564	P759362-091	P750382-091	P114492-091		P764321-469
	S/S	P764315-581	P764315-582	P764315-583	P764315-583	P764315-584		P764321-469	
20AS Solid State* (3) 6.5 KW	C/S	P764319-903	P764319-903	P764319-903	P764319-903	P764319-916	P764329-282		P764317-713 (Set of 3)
	S/S	P764321-272	P764321-272	P764321-272	P764321-272	P764321-398			
30AS Solid State* (3) 9.75 KW	C/S	P764322-460	P764322-460	P764322-460	P764322-460	P764322-882	P764329-283		
	S/S	P764322-394	P764322-394	P764322-394	P764322-394	P764322-395			
LB 20 Electric (1) 20 KW Indirect (Stm to Stm)	C/S	P078826-091	P075425-091	P764317-564	P759362-091	P750382-091	P114492-091		P764321-472 (Set of 3)
	S/S	P764315-581	P764315-582	P764315-583	P764315-583	P764315-584		P764321-344	
LB 30 Electric (1) 30 KW Indirect (Stm to Stm)	C/S							P764324-729	
	S/S	P758665-091	P753750-091	P759940-001	P761959-001	P753995-091	P114650-091		
	C/S	P764321-563	P764321-564	P764321-565	P764321-565	P764321-566			
	S/S							P764321-344	
	S/S							P764324-729	

\*Non-Solid State = Gray Cabinet  
Solid State = Pump & Motor on Top of Generator cabinet.

C/S = Carbon Steel  
S/S = Stainless Steel

NOTE: 1. All heaters include heater gasket.  
2. All stand alone LB Series generators have gray

# ELECTRO-STEAM STEAM GENERATOR HEATING ELEMENTS

Size/Type	Mat'l	Voltage						Indirect (Stm to Stm)	Heater Gasket
		208	240	380	415	480	575/600		
LB40	Electric	P078826-091	P075425-091	P764317-564	P759362-091	P750382-091	P114492-091		
	(2) 20 KW	P764315-581	P764315-582	P764315-583	P764315-583	P764315-584			
	Indirect						P764321-344		
	(Stm to Stm)						P764324-729		
LB60	Electric	P758665-091	P753750-091	P759940-001	P761959-001	P753995-091	P114650-091		
	(2) 30 KW	P764321-563	P764321-564	P764321-565	P764321-565	P764321-566			
	Indirect						P764321-344		P764321-472 (Set of 3)
	(Stm to Stm)						P764324-729		
LB80	Electric	P759377-091	P753854-091	P764316-029	P764316-029	P753994-091	P115851-000		
	(3) 25 KW	P764319-441	P764319-442	P764316-031	P764319-443	P764319-444			
	Indirect						P764321-344		
	(Stm to Stm)						P764324-729		
LB100	Electric	P759377-091	P753854-091	P764316-029	P764316-029	P753994-091	P115851-000		
	(4) 25 KW	P764319-441	P764319-442	P764316-031	P764319-443	P764319-444			
	Indirect						P764321-344		
	(Stm to Stm)						P764324-729		

C/S = Carbon Steel  
S/S = Stainless Steel

NOTE: 1. All heaters include heater gasket.  
2. All stand alone LB Series generators have gray cabinet.



# **ELECTRO-STEAM STEAM GENERATOR HEATING ELEMENTS**

Size/Type	Mat'l	Voltage						Indirect (Stm to Stm)	Heater Gasket
		208	240	380	415	480	575/600		
LB 120	Electric	P758665-091	P753750-091	P759940-001	P761959-001	P753995-091	P114650-091		P764321-472 (Set of 3)
	(4) 30 KW	P764321-563	P764321-564	P764321-565	P764321-565	P764321-566			
	Indirect							P764321-344	
	(Stm to Stm)							P764324-729	
LB 150	Electric	P759377-091	P753854-091	P764316-029	P764316-029	P753994-091	P115851-000		
	(6) 25 KW	P764319-441	P764319-442	P764316-031	P764319-443	P764319-444			
	Indirect							P764321-344	
	(Stm to Stm)							P764324-729	
LB 180	Electric	P758665-091	P753750-091	P759940-001	P761959-001	P753995-091	P114650-091		P764321-472 (Set of 3)
	(6) 30 KW	P764321-563	P764321-564	P764321-565	P764321-565	P764321-566			
	Indirect							P764321-344	
	(Stm to Stm)							P764324-729	

C/S = Carbon Steel  
S/S = Stainless Steel

NOTE: 1. All heaters include heater gasket.  
2. All stand alone LB Series generators have gray cabinet.

# CHROMALOX STEAM GENERATOR HEATING ELEMENTS

Size/Type	Mat'l	Voltage							Heater Gasket
		208	240	380	415	480	575	600	
CHS 20AS	C/S	P764325-727	P764325-727	P764326-982	P764326-982	P764325-726			P764323-527
CHS 30AS	C/S	P764323-549	P764323-549	P764323-639	P764323-639	P764323-550	P764327-911	P764328-689	
	S/S	P764324-764	P764324-764	P764327-659	P764327-659	P764324-765			
CHS 30AS Century	C/S	P764328-020	P764323-549	P764323-639	P764323-639	P764323-550	P764327-911	P764328-689	P764323-527
After date code 304LP	S/S	P764328-094	P764324-764	P764327-659	P764327-659	P764324-765			
CES 20 (1) 20 KW	C/S	P764323-413	P764323-415			P764323-417			
(CH04)	S/S	P764326-594							
CES 30 (1) 30 KW	C/S	P764323-418	P764323-419	P764325-448		P764323-420			
(CH05)	S/S	P764326-598	P764326-753	P764327-659		P764326-754			
CES 40 (2) 20 KW	C/S	P764324-886	P764323-424	P764324-887		P764323-425			P764323-594
(CH06)	S/S	P764326-596	P764325-671						
CES 60 (2) 30 KW	C/S	P764323-418	P764323-419	P764325-448		P764323-420			
(CH07)	S/S	P764326-598	P764326-753	P764327-659		P764326-754			
CES 75 (2) 37.5 KW	C/S	P764323-421	P764323-422	P764326-083	P764328-077	P764323-423			
(CH08)	S/S		P764326-300			P764326-364			
CES 100 (3) 34 KW	C/S	P764323-655	P764323-657	P764324-731	P764326-875	P764323-691		P764329-304	
(CH10)	S/S		P764326-597			P764323-999			
*575V = P764327-911 600V = P764328-689		C/S = Carbon Steel S/S = Stainless Steel			NOTE: 1. All heaters include heater gasket. 2. All stand alone CH Series generators have beige cabinet.				

NOTE: 1. All heaters include heater gasket.  
2. All stand alone CH Series generators have beige cabinet.

C/S = Carbon Steel  
S/S = Stainless Steel

\*575V = P764327-911  
600V = P764328-689

# CHROMALOX STEAM GENERATOR HEATING ELEMENTS

Size/Type	Mat'l	Voltage						Heater Gasket
		208	240	380	415	480	575/600	
CES 120 (3) 40 KW (CH12)	C/S	P764324-889	P764325-065	P764326-866	P764326-867	P764323-715		P764323-594
	S/S	P764329-363				P764326-595		
CES 150 (4) 37.5 KW (CH13)	C/S	P764323-421	P764323-422	P764326-083	P764328-077	P764323-423		
	S/S		P764326-300			P764326-364		
CES 180 (4) 45 KW (CH14)	C/S	P764326-541	P764323-426	P764328-076	P764328-075	P764323-427		
	S/S					P764327-774		

C/S = Carbon Steel  
S/S = Stainless Steel

NOTE: 1. All heaters include heater gasket.

2. All stand alone CH Series generators have beige cabinet.

## SERVICE SUPPORT RESOURCES ERIE CONTACTS

Function	Phone Ext.	Fax Number
<u>Contract Department:</u> V. Panetta	2641	(814) 870-8841
<u>Service Support:</u> M. Reske	2509	800-270-7624
<u>Parts Ordering:</u> Hunt Group	3070	800-270-7624
<u>Inventory Control:</u> K. Weibel	2591	800-482-0056
<u>Field Service Engineering:</u> Hunt Group Toll Free Number	3053 1-800-825-0852	(814) 870-8201
<u>Return Goods/Calibration:</u> J. Gennuso	37954	(440) 392-7954
<u>Specialty Service:</u> M. Major	3503	(814) 870-8201
<u>Specialty Support/Pricing:</u> R. Torok	2507	(814) 870-8453
<u>Service Training:</u> M. Varhola	8163	(814) 870-8164

## KEY ACCOUNTS

Account - A/R \_\_\_\_\_

Contact \_\_\_\_\_

Phone No. \_\_\_\_\_

Account - A/R \_\_\_\_\_

Contact \_\_\_\_\_

Phone No. \_\_\_\_\_

Account - A/R \_\_\_\_\_

Contact \_\_\_\_\_

Phone No. \_\_\_\_\_

Account - A/R \_\_\_\_\_

Contact \_\_\_\_\_

Phone No. \_\_\_\_\_

Account - A/R \_\_\_\_\_

Contact \_\_\_\_\_

Phone No. \_\_\_\_\_

Account - A/R \_\_\_\_\_

Contact \_\_\_\_\_

Phone No. \_\_\_\_\_

## KEY ACCOUNTS

Account - A/R \_\_\_\_\_

Contact \_\_\_\_\_

Phone No. \_\_\_\_\_

Account - A/R \_\_\_\_\_

Contact \_\_\_\_\_

Phone No. \_\_\_\_\_

Account - A/R \_\_\_\_\_

Contact \_\_\_\_\_

Phone No. \_\_\_\_\_

Account - A/R \_\_\_\_\_

Contact \_\_\_\_\_

Phone No. \_\_\_\_\_

Account - A/R \_\_\_\_\_

Contact \_\_\_\_\_

Phone No. \_\_\_\_\_

Account - A/R \_\_\_\_\_

Contact \_\_\_\_\_

Phone No. \_\_\_\_\_

## NOTES

[illegible]

[illegible]